

3 IN. X 3 IN. DENTAL MIXING PAD

REV022522

Component of the Dental Emergency Response Kit

A Henry Schein Brand extra-heavy-duty poly pads features a high-gloss finish that is great for mixing cements and composites. With a non-skid-foam base. Poly mixing pads ensure a clean, smooth mix without lint or paper shreds.

Specifications:

- 3 in. X 3 in. extra-heavy-duty
- High gloss finish
- Non-Skid foam base
- Great for dental cements, composites, and liners

Dimensions:

- L 3 in. x W 3 in.
- Weight: 0.11 lbs

3 in. x 3 in. Dental Mixing Pad

ITEM #	QUANTITY	NSN #
ZZ-1421	1 Each	--



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®