

Products available on DAPA

D-BCRK *Department of Defense Bleeding Control Response Kits - Individual (Nylon or Vacuum-Sealed)*

Control the bleed with these easy-to-use tools!

Life threatening bleeding can occur from everyday emergencies such as vehicle crashes, incidents involving heavy equipment, accidental injuries from natural disasters like tornadoes, and inflicted injuries from man-made disasters such as shootings and bombings. D-BCRK Individual Bleeding Control Kits by North American Rescue are designed to provide essential equipment that empowers DOD Forces and Federal Employee responders to take immediate action in stopping life threatening bleeding.

Contained in a ruggedized, nylon carry bag or vacuum-sealed pouch, this life saving kit provides intuitive and easy-to-use tools that are proven to help to save lives such as the fastest, safest and most effective prehospital tourniquet in the world, the C-A-T®, quality gauze, shears, and pressure dressings/bandages. D-BCRK Individual Bleeding Control Kits are available in Basic, Intermediate, and Advanced configurations.

Nylon Bag Option:

- Easy to restock in cases of use
- Convenient to carry and easy to open in stressful situations
- Opens flat for ease of use and immediate identification of enclosed devices
- Easily recognizable red nylon bag with ample hook & loop field for attaching an ID patch

Vacuum Sealed Option:

- Cost-effective alternative to Nylon
- Compact and rugged Vacuum-Sealed pouch
- Large Identifying Card Vacuum-Sealed inside the pouch for ease of identification

Specifications:

- Nylon Bag: H 7 in. x W 5 in. x D 2.75 in.
- Vacuum Sealed: H 7 in. x W 5 in. x D 2.75 in.

Every **D-BCRK INDIVIDUAL BLEEDING CONTROL KIT** contains:



QTY 1 C-A-T® Tourniquet
 QTY 1 ETD™ 6 in. Dressing
 QTY 1 Gauze S-Rolled
 PAIR 2 Nitrile Gloves, Lg
 QTY 1 Trauma Shears, Lg

PLUS these Individual Components based upon configuration of choice:



BASIC
 QTY 1 Gauze S-Rolled

INTERMEDIATE
 QTY 1 Gauze S-Rolled
 QTY 1 HyFin® Vent, Twin Pk

ADVANCED
 QTY 1 Combat Gauze® Z-fold
 QTY 1 HyFin® Vent, Twin Pk



Vacuum Packed or in a Nylon bag with hook & loop field for attaching an ID patch

REV091720

D-BCRK Bleeding Control Response Kits - Individual

LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	15 oz	85-2057
Intermediate	Red Nylon	16.5 oz	85-2058
Advanced	Red Nylon	17.4 oz	85-2059
Basic	Vacuum Sealed	11.7 oz	85-1749
Intermediate	Vacuum Sealed	13.25 oz	85-2066
Advanced	Vacuum Sealed	14 oz	85-2069



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

D-BCRK Individual Bleeding Control Response Kits also available in 8-packs and Wall Stations (see other side for more details)

Products available on DAPA

D-BCRK Department of Defense Bleeding Control Response Kits - 8-Packs & Stations (Nylon or Vacuum-Sealed)

REV091720

Stop life-threatening bleeding in the first critical minutes...

D-BCRK (Department of Defense Bleeding Control Response Kits) 8-Packs and Wall Stations from North American Rescue are designed to provide essential equipment that empowers DOD Forces and Federal Employee responders to take immediate action in stopping life threatening bleeding. Each 8-Pack/Wall Station includes eight (8) Individual Bleeding Control Kits (Vacuum-packed or with Red Nylon bag /PVC Bleeding Control patch) that are contained inside of an easy-to-open nylon carrying case. The wall stations also include a durable, clear polycarbonate cabinet designed to safely store the D-BCRK Bleeding Control 8- Pack on a wall in public areas.

D-BCRK Bleeding Control Response Kits | 8-Packs and Wallstations are available in Basic, Intermediate, and Advanced configurations.

D-BCRK 8-PACK Packaging Options:



See other side for contents of each Individual Bleeding Control Kit

D-BCRK CLEAR WALL STATIONS also include:



Polycarbonate Wall Case:

- Durable, clear polycarbonate cabinet with rounded corners provides both physical and environmental protection for medical equipment
- Includes a wire shelf to hold the 8-pack & an external eyelet that allows for a secondary safety seal to easily identify if the station has been opened or tampered with
- Clear design allows visibility at any angle
- External mounting tabs and hardware enable the case to be easily installed on any wall
- Alarm available upon special request

8-Packs:

- Eight (8) Nylon or Vacuum-sealed Individual Bleeding Control Kits contained inside of an easy-to-open nylon carrying case
- Available in Basic, Intermediate, and Advanced configurations

Specifications:

Clear Wall Case:

- H 18 in. x W 15.5 in. x D 10 in.

Each 8-Pack:

- H 14 in. x W 11 in. x D 8 in.

Individual Bleeding Control Kit:

- H 7 in. x W 5 in. x D 2.5 in.

D-BCRK Bleeding Control Response Kits - 8-Pack

LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	8.5 lb	85-2072
Intermediate	Red Nylon	9.3 lb	85-2073
Advanced	Red Nylon	9.7 lb	85-2074
Basic	Vacuum Sealed	8.5 lb	85-1750
Intermediate	Vacuum Sealed	9.3 lb	85-2067
Advanced	Vacuum Sealed	9.7 lb	85-2070

D-BCRK Bleeding Control Response Kits - Station

LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	23.5 lb	85-2075
Intermediate	Red Nylon	24.3 lb	85-2076
Advanced	Red Nylon	24.7 lb	85-2077
Basic	Vacuum Sealed	23.5 lb	85-1751
Intermediate	Vacuum Sealed	24.3 lb	85-2068
Advanced	Vacuum Sealed	24.7 lb	85-2071



#ZZ-0822
Clear Polycarbonate Wall Case Only also available for separate purchase



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277