

# MEDIC TRAUMA PACK

REV020819

*Light Framed Platform and Removable Pull Bags for all your Trauma Management Essentials*



## **Our *Most Comprehensive* Trauma Management Platform for Tactical Healthcare Professionals in the Line of Duty**

The Medic Trauma Pack is a light framed platform with silent, dirt and sand-resistant zippers that close over a large main compartment with modular, removable pull bags for all your trauma management essentials.

The adjustable shoulder and sternum straps make for a comfortable fit. Cinch straps are included for the sides to maintain load integrity during stress maneuvers. Grommets are provided at the base of the pack for drainage.

Includes 9 removable, color-coded pull bags and 2 removable side compartments to provide modular organization for all your medical equipment. With over 1,000 cubic inches of gear capacity this is our most comprehensive kit/pack available.

### **Specifications:**

- Case Closed: H 19 in. x W 19 in. x D 10 in.
- Weight: 18 lbs 1 oz

### **Medic Trauma Packs**

Color	Item#	NSN#
BLK	80-0070	6545-01-572-7251
ODG	80-0071	6545-01-594-5617*
COY	80-0072	6545-01-588-1862
DUC	80-0096	

### **Available Colors:**



\*NATO NSN#



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
— WITH A —  
**MISSION®**

# MEDIC TRAUMA PACK

REV020819

## Features:

- Over 1000 cubic inches of gear capacity
- Contoured, padded shoulder straps to distribute load over a larger area for comfortable carry
- Flap allows for Hydration Tube to extend out of pack
- Constructed of durable 500 Denier Nylon
- Sand & water resistant zipper with silent cord pulls
- External Cinch Straps to tighten pack and decrease load shift
- MOLLE Style Attachment System Compatible

## Kit Contents:

- 1 x Bag (CCRK® - Medic Trauma)
- 5 x pair Black Talon® Nitrile Trauma Gloves
- 2 x NAR PPE Kit™ (Personal Protection Equipment Kit)
- 2 x Nasopharyngeal Airway 28F with Lubricant
- 1 x King LT-D™ Supraglottic Airway
- 1 x Tactical Suction Device
- 1 x Cyclone® New & Improved BVM
- 2 x Petrolatum Gauze (3 in. x 18 in.)
- 2 x ARS® Needle Decompression Kit (14 G x 3.25 in.)
- 2 x C-A-T® (Combat Application Tourniquet®)
- 6 x ETD™ 6 in. Emergency Trauma Dressing
- 1 x NAR Sharps Container
- 6 x Polycarbonate Eye Shield
- 4 x S-Rolled Gauze™ (4.5 in. x 4.1 yd)
- 2 x ETD™ Abdominal Emergency Trauma Dressing
- 4 x NAR Saline Lock Kit
- 1 x BOA® IV Constricting Band
- 2 x Surgical Tape (2 in.)
- 1 x Trauma Shears (7.25 in.)
- 4 x SAM® Splint II
- 1 x TTS™ Tactical Traction Splint
- 1 x Heat Reflective Shell
- 1 x Armadillo™ Medication Storage Case
- 1 x Combat Casualty Reference Card
- 4 x Triage Card (NAR T2 Tag®)
- 1 x NAR Headlamp

## MEDIC TRAUMA PACK:



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**