

NAR-4™ CHEST POUCH

REV011419



Streamlined to provide Tactical Healthcare Professionals rapid access...

The auxiliary NAR-4™ companion is streamlined to provide Tactical Healthcare Professionals with rapid access to critical basic & advanced equipment on their chest.


The container is constructed of 500D CORDURA® Solution-Dyed Nylon with Near Infrared Signature Reduction, Ghillie-Tex IR signature reduction hardware and a silent, self repairing nylon zipper with 550 pull cords.

The kit is designed to be secured to the front of the body armor using a MOLLE-style attachment system.

Specifications:

- Cased Closed: L 7.75 in. x W 9.75 in. x 4.5 in.
- Weight: 2 lb 10.3 oz

Available Colors:



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
— WITH A —
MISSION®


NAR-4™ CHEST POUCH

REV011419

Kit Contents:

- 1 x Bag (CCRK® - NAR-4™ Chest Pouch)
- 3 x Black Talon® Nitrile Glove Kit, Lg.
- 2 x Combat Casualty Card
- 2 x Nasopharyngeal Airway 28F with Lubricant
- 1 x HyFin® Vent Chest Seal Twin Pack
- 2 x ARS® Needle Decompression Kit (14 G x 3.25 in.)
- 2 x C-A-T® (Combat Application Tourniquet®)
- 2 x ETD™ 6 in. Emergency Trauma Dressing
- 2 x S-Rolled Gauze™ (4.5 in. x 4.1 yd)
- 2 x Choice of Hemostatic (Combat Gauze or Chito Gauze PRO)
- 1 x Scissor Leash™
- 1 x Trauma Shears (7.25 in.)
- 1 x Permanent Marker
- 1 x Tactical Crickit®
- 2 x Polycarbonate Eye Shield

Color	Hemostatic	Item#	NSN#
BLK	ChitoGauze	80-0283	
COY	ChitoGauze	80-0284	
DUC	ChitoGauze	80-0285	
MTC	ChitoGauze	80-0287	
ODG	ChitoGauze	80-0288	
BLK	Combat Gauze	80-0173	
COY	Combat Gauze	80-0175	
DUC	Combat Gauze	80-0176	
MTC	Combat Gauze	80-0247	
ODG	Combat Gauze	80-0289	



Features:

- Configurable based on individual skill set and mission requirements
- Rapid access and organization is enhanced by elastic storage loops that keep items in place until needed
- 500D CORDURA® Solution-Dyed Nylon with built in Near Infrared Signature Reduction
- Sand and water resistant zippers with silent pulls
- PALS/MOLLE-style attachment system
- Ambidextrous Trauma Shear configuration



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®