

USMC HYFIN® COMBO PACK *Chest Seals*

REV112519

The USMC HyFin® Chest Seal Combo Pack is a versatile kit designed for the treatment of both entry/exit or multiple penetrating injuries to the chest. Each combo pack contains one (1) HyFin® Occlusive Chest Seal and one (1) HyFin® Vent Chest Seal with 3-channel pressure relief vents all packaged in a convenient clear resealable storage pouch. Both seals feature enhanced adhesive technology designed for a superior seal in the most adverse conditions, including sweaty or hairy casualties. In addition, their transparent backing facilitates placement directly over the wound and easily conforms to the casualty's chest. The large, easily-visible red elliptical tabs are designed to allow for single step peel-and-apply application and to facilitate the burping of a wound if necessary.

Features:

- Versatile combo pack designed for the treatment of both entry/exit or multiple penetrating injuries to the chest
- Featuring NAR's patented HyFin® Chest Seals with enhanced adhesive technology designed for a superior seal in the most adverse conditions, including sweaty or hairy casualties
- Includes one HyFin® Occlusive & one HyFin® Vented Chest Seal packaged in a clear resealable pouch
- Large, easily-visible red elliptical tabs designed to allow for single step peel-and-apply application and the burping of a wound if necessary

Dimensions:

- Packaged: L 9 in. x W 5 in. x D 0.2 in.
- Each Chest Seal (when deployed): H 6 in. x W 6 in. x D 0.1 in.
- Total Weight: 2.9 oz
- US Patent 7,504,549 & Patent(s) Pending

Combo Pack Contents:

- 1 x HyFin® Chest Seal (Occlusive)
- 1 x HyFin® Vent Chest Seal (with 3-channel pressure relief vents)
- 1 x Resealable clear pouch



USMC HyFin® Chest Seal Combo Pack

| Item # | NSN# |
|---------|------------------|
| 85-1045 | 6510-01-654-6787 |



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**