

JETT Junctional Emergency Treatment Tool

REV120318

Simple and Safe Solution for Occluding Blood Flow to both Lower Extremities...

The Junctional Emergency Treatment Tool (JETT™) includes two pressure devices in a single tool, allowing for simultaneous occlusion of blood flow to both lower limbs (bilateral) with an easy to apply, pre-assembled, ready-to-use belt solution. The unique design of the JETT™ allows for patient transport while the device is in place and includes a lanyard and toggle solution to ensure that the Windlass Axle does not loosen due to vibration or bumps during transport.

In cases where there are bilateral injuries, the JETT™ offers a simple, safe solution for occluding blood flow to both lower extremities. This expanded capability is built into every JETT™, meaning that with standardized application training the device will be pre-positioned for rapid bilateral deployment whether applying for a unilateral or bilateral injury.

High extremity wounds in the thigh/groin region where tourniquets cannot be placed or are less effective require equipment such as the JETT™ to stop the bleeding and potentially save the life of a casualty that would otherwise bleed out in minutes. Penetrating trauma, blast injuries and amputation are examples of injuries that may require treatment with this device.

JETT - Junctional Emergency Treatment Tool

Item #	NSN#
30-0088	6515-01-616-5841

Stop bleeding and potentially save the life of a casualty that would otherwise bleed out in minutes...



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®

JETT *Junctional Emergency Treatment Tool*

REV120318

Features:

- Compact, lightweight & ruggedized for use in the tactical environment
- Pre-assembled, Ready-to-use & Easy-to-apply at or near the point of injury
- Simultaneous Occlusion of blood flow to both lower limbs that does not impede Respiration
- 2 Compression Pads, pre-positioned and individually adjustable, allow for individual or bilateral activation, based on the patient's injuries
- Circumferential Compression Buttocks Pad safely secures the device to the patient
- Lanyard & Toggle Locking Mechanism ensures that vibrations or bumps to the windlass will not reduce the pressure during patient movement
- Writeable Area to document the time of application
- Easy-to-understand, Step-by step Instructions included
- Sleeves help pre-position the pads & aid in smooth application

Dimensions:

- Packaged: H 5.5 in. x W 7 in. x D 3.75 in.
- Overall Device Length: 56.75 in. Fits waist up to 50 in.
- Weight: 1 lb 8 oz

Nylon pouch designed for carrying the JETT® on the outside of an aid bag is also available (#80-0355 Black #80-0356 Coyote)



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**