

# REFLEX VEHICLE MOUNT

REV101323

The ReFlex™ Vehicle Mount is a companion accessory to the ReFlex™ IFAK System, specifically designed to mount the ReFlex™ Med Roll to a vehicle head rest. It features an adjustable, nylon-and-elastic cinch strap closure that secures with a 1" side-release buckle, enabling this to fit the Med Roll regardless of how full it is packed. This mount attaches with two 1", adjustable nylon webbing straps that secure by threading through a tri-glide. The Med Roll can be removed quickly by pulling to the side, without the need to unbuckle the retention strap. This mount also includes a 3" elastic loop for the addition of a tourniquet.

## Specifications:

- Attaches ReFlex™ Med Roll to the back of a vehicle headrest
- Universal design fits virtually all vehicles with detachable headrests
- Elastic cinch strap allows instant retrieval of the Med Roll from either side of the vehicle
- Built-in elastic tourniquet loop
- All-nylon construction

## Dimensions:

- L 5.5 in. x W .125 in. x H 6.0 in.
- Weight: 0.21

## Reflex Vehicle Mount

ITEM #	COLOR	NSN #
ZZ-1369	Coyote	
ZZ-1370	Gray	



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**