

FIRST LIGHT TOMAHAWK TC3 NV *Night Vision Tactical Lights*

REV092322

Specifically Engineered to assist the Combat Medic...

In response to the growing, and to this point, unmet need from tactical medical professionals, NAR is proud to offer the latest innovation from First Light USA® - the industry leader in combat-specific lighting tools. The Tomahawk® Night Vision (NV) Series of hands-free tactical lights are specially designed to provide maximum versatility and enhanced illumination at the point of injury.

The hands-free design adds a unique set of features specifically engineered to assist the Combat Medic— from the superior chest mounting capability and full 180° adjustability of the light source to the proprietary mix of red and green LEDs that “highlight” blood for immediate wound identification in multi-fluid situations.

This not only provides Medics with a quick tool to accurately assess a casualty in low light/night time conditions, it allows them to maintain light discipline as well as a defensible position.

Features:

- Right angle tactical light with rotating head and 3 button control panel for one-handed thumb actuation
- Provides White, IR, and Red+Green LED mix for highlighting blood during casualty treatment
- Featuring Integrated high intensity IR signaling beacon
- Includes momentary and three intensity levels of constant-on (all output colors)
- Provides 500 lumens of maximum white light, 2 hour Runtime on high output
- Can be used handheld or hands-free with MOLLE TRS Blade Mount
- Lock-out mode prevents accidental illumination

Specifications:

- Output Wavelengths: White, IR, & Red+Green (Blood Highlighting), IR Signaling Beacon
- Maximum White Light Output: 425 lumens
- Max Output Runtime: 2 hours (ANSI Tested)
- Power Source: CR123 Batteries (2)
- Material: Aerospace-Grade Aluminum

Dimensions:

- With mount: H 3.4 in. x L 2.4 in. x W 1.5 in.
- Weight (with mount & batteries): 8 oz

In the Box:

- 1 x Right angle tactical light with rotating head and 3 button control panel
- 2 x Retention System Finger loops
- 1 x MOLLE Tactical Retention System (TRS) Blade Mount
- 2 x CR123 batteries
- 1 x User Manual
- 1 x Task Light (UV or Blue)



First Light Tomahawk TC3 Night Vision Tactical Light

COLOR	TASK LIGHT	ITEM#
Coyote	UV	70-0015
OD Green	UV	70-0112
Coyote	Blue	70-0111
OD Green	Blue	70-0110



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH

MISSION®