

# NAR-5 SAR AID KIT™

REV110118

**The standard casualty response bag which meets the requirements for the SAR Level A Medical Bag**

The NAR 5 Aid Kit™ from North American Rescue is designed to meet the rigorous requirements of the U. S. Navy Search and Rescue (SAR) Swimmers. This ergo dynamic, low profile & durable casualty response bag is constructed of Berry Amendment 500D CORDURA® solution dyed nylon with GhillieTEX™ hardware for near IR signature reduction and sand/water resistant, self repairing zippers with silent 550 cord pulls and is **designed to be hung in aircraft or other transport vehicle by top bag handle**. The NAR 5 also features 6 hook-and-loop internal pouches in various sizes which can be self labeled for the specific mission set; 4 smaller pouches for hemorrhage control products and 2 larger ones to accommodate airway & intravenous products.

Pouches are designed for quick access to medical equipment necessary to respond to specific life threatening medical emergencies. Pouches attach via hook-and-loop to allow the first responder to remove the specific pouch needed to provide on-scene medical care to the patient.

Slim-line, low-profile bag is securable using anti-pilferage seals to protect all medical components inside the bag. Kit can be carried as a backpack using shoulder straps that stow away until needed.

## Features:

- Meets the requirements for the SAR Level A Medical Bag
- Ergo dynamic, low profile & durable casualty response bag constructed of Berry Amendment 500D CORDURA® solution dyed nylon
- GhillieTEX™ hardware for near IR signature reduction and sand/water resistant, self repairing zippers with silent 550 cord pulls
- Featuring 6 hook-and-loop internal pouches in various sizes which can be self labeled for the specific mission set; 4 smaller pouches for hemorrhage control products and 2 larger ones to accommodate airway & intravenous products

## Specifications:

- Dimensions: H 24 in. x W 15 in. x D 7 in.
- Kit Weight: 14 lb 10 oz

*Designed to be hung in aircraft or other transport vehicle by top bag handle*

COY

## NAR 5 Aid Kit

Options	Item#	NSN#
Aid Kit	80-0416	
Bag Only	80-0399	8465-01-669-5796*

\*NATO NSN



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**

# NAR-5 SAR AID KIT™

REV110118



## Kit Contents:

- 1 x Bag (CCRK® - NAR-5 SAR Aid)
- 4 x C-A-T® (Combat Application Tourniquet®)
- 4 x NAR Compressed Gauze
- 4 x ETD™ 4 in. Emergency Trauma Dressing
- 1 x ETD™ Abdominal Emergency Trauma Dressing
- 4 x HyFin® Vent Twin Pack Chest Seal
- 1 x Cyclone® BVM
- 1 x CPR Mask
- 4 x Triangular Bandage
- 4 x 3 in. Elastic Wrap Bandage
- 2 x SAM® Splint II
- 2 x Cervical Collar
- 1 x Roll Surgical Tape (1 in.)
- 1 x Trauma Shears (7.25 in.)
- 1 x Sterile Burn Dressing
- 2 x 1000cc IV Solution
- 4 x NAR PPE Kit™ (Personal Protection Equipment Kit)
- 12 x Flexible Fabric Bandages
- 1 x Survival Blanket
- 1 x Mini Maglight with batteries



*Designed to meet the rigorous requirements of the U. S. Navy Search and Rescue (SAR) Swimmers*



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**