

# SIMULATOR BLOOD FILL TANK

REV011624

## ***Washable and Reusable***

The Simulator Blood Fill Tank is a 2.5 gallon manually pressurized reservoir designed to pump the simulated blood concentrate into North American Rescue Simulation® training simulators.

## **Specifications:**

- To be used with all TOMManikins, TAMkin, and K9 simulators

## **Dimensions:**

- Packaged - L 36 in. x W 18 in. x H 18 in.
- Weight: 4 lbs

## ***Spare Connectors***

ITEM #	NSN #
93-0183	----



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**