

# T.A.L.O.N.™ MANUAL NEEDLE SET (BOX OF 6)

BY EZ-IO®

REV061526



## Manual I/O Insertion...

Manual I/O insertion at all 7 extremity sites, including the sternum. Includes 15 gauge, 304 stainless steel Needle/Catheter, Sharps Protector, EZ-Connect® Extension Set, I/O Wristband, Quick Reference Card & Sternal Locator. Sterile.

\*\*\* EZ-IO® Products are restricted to Federal Purchases only. \*\*\*

\*\*\* Please contact Customer Service for details. \*\*\*

## Features:

- Sold in Box of 6
- The EZ-IO® T.A.L.O.N.™ Manual Needle Set allows for manual I/O insertion at all 7 extremity sites (sternum, bilateral proximal humerus, bilateral proximal tibia, and bilateral distal tibia)
- Includes 15 gauge, 304 stainless steel Needle/Catheter, Sharps Protector, EZ-Connect® Extension Set, I/O Wristband, Quick Reference Card & Sternal Locator

## Specifications:

- Packaged: L 10 in. x W 8.5 in
- Weight: 6.5 oz

## T.A.L.O.N.™ Manual Needle Set

ITEM #	NSN #
30-0132	6515-01-642-7919



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**