

OCHO MESH AID POUCH

GEAR ORGANIZERS



Additional Storage and Rapid Access...

The Ocho Mesh Aid Pouch from North American Rescue is a customizable, modular solution for outfitting your low visibility bag with the right medical gear for any mission.

With its hook-and-loop backing, zippered mesh pouch, and red tip pull tabs, the pouch keeps equipment organized and securely intact until needed.

Designed to provide rapid access & additional organizational storage for all your critical gear, the Ocho family of quick detach modules are the convenient solution to keep everything you need right at your fingertips.

Features:

- Mesh pouch features a detachable label for easy identification of hemorrhage control equipment (sold separately)
- Provides quick detach capability using the red quick-pull tab
- Constructed of 500D CORDURA® Berry Amendment Nylon, solution dyed for near IR signature reduction
- Backed with milspec hook and compatible with any storage solution with a loop panel attached
- *Medical equipment sold separately*

Specifications:

- L 7 in. x W 4.5 in. x D 3 in.
- Weight: 1.87 oz



Color	Item#	NSN#
BLK	80-0259	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS WITH A
MISSION®