

# CASUALTY IMMOBILIZATION SYSTEM

*Easy to Use Patient Strapping System to Help Secure Casualties to Spine Boards*

The Casualty Immobilization System is constructed of heavy duty polypropylene webbing for fast and effective immobilization. Its unique design comes with six points of attachment and secures easily to the shoulder and foot area of the spine board.

## CASUALTY IMMOBILIZATION SYSTEM

ITEM #	SEL #	NSN #
50-0031	09ME-04-SPIN	6530-01-586-2411

## FEATURES/SPECIFICATIONS

- Total backboard strap system for fast, firm total body immobilization
- Two "suspender" type over-the-shoulder straps
- No choking "V" by patients neck
- 4 adjustable torso/limb straps, with adjustable center strap
- Black polypropylene webbing with hook and loop closures
- *NOTE: All accessories are sold separately. Spine Boards are not included.*



**Six points of attachment and secures easily to the shoulder and foot area of the spine board.**

## DIMENSIONS

- L 16.5 in. x W 4 in. x 2.75 in.
- Weight: 1 lb 0.5 oz



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277