

NAR CASUALTY IMMOBILIZATION System

TACEVAC

Easy to Use Patient Strapping System to Help Secure Casualties to Spine Boards...

The NAR Casualty Immobilization System is constructed of heavy duty polypropylene webbing for fast and effective immobilization. Its unique design comes with six points of attachment and secures easily to the shoulder and foot area of the spine board.

Six points of attachment and secures easily to the shoulder and foot area of the spine board.

Features:

- Total backboard strap system for fast, firm total body immobilization
- Two "suspender" type over-the-shoulder straps
- No choking "V" by patients neck
- 4 adjustable torso/limb straps, with adjustable center strap
- Black polypropylene webbing with hook and loop closures

Specifications:

- L 16.5 in. x W 4 in. x 2.75 in.
- Weight: 1 lb 0.5 oz

SEL#:
SEL#: 09ME-04-SPIN

NOTE: All accessories are sold separately.
Spine Boards are not included.

NAR Casualty Immobilization System

Item#	NSN#	CE
50-0031	6530-01-586-2411	

Total backboard strap system for fast, firm total body immobilization.

PRODUCTS WITH A
MISSION[®]



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277