

# BIODEGRADABLE HRP *Human Remains Pouch*

REV070523

The Biodegradable Human Remains Pouch (HRP) from North American Rescue offers a sturdy solution for containing and moving remains. The HRP is lightweight and comes packaged in a rugged vacuum sealed bag that minimizes the storage space requirements. The pouch features 8 large carry handles that enable multiple people to work together when moving. The HRP meets OSHA REG 3130 for the retention of bodily fluids and the 1.5 inch webbing support straps offer a weight carry capacity of up to 350 lbs.

## Features:

- Large curved zipper for easy access
- Light weight design
- Packaged in vacuum seal bag for minimum cube space
- 8 large carry handles for multiple person carry
- 1.5" webbing support straps with integrated carry handles
- Designed to contain bodily fluids per OSHA REG 3130
- Weight carry capacity up to 350 lbs

## Specifications:

- Weight: 3.4 lb
- Product Opened: H 94 in. x W 36 in. x D 0.25 in.
- Packaged: H 14 in. x W 10 in. x D 2.5 in

### **Biodegradable HRP - Human Remains Pouch**

ITEM#	NSN:
40-0003	6530-01-492-3169



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH  **MISSION®**