


ADC Stethoscope

Full-range Stethoscope for Tactical Health Care Professionals...

This full-range Stethoscope is for Tactical Health Care Professionals who demand the best in auscultation and comfort.

This is the most comfortable stethoscope manufactured today. Superior acoustics mean accurate auscultation of heart, lung, and blood pressure sounds.



*Best full-range
stethoscope for
auscultation
and comfort.*

STETHOSCOPE BY ADC

ITEM #	NSN
20-0060	----

FEATURES/SPECIFICATIONS

- Full-range stethoscope for Tactical Health Care Professionals who demand the best in auscultation and comfort
- Adult diaphragm disc with non-chill rim
- Double-head chest-piece assembly

- Hard-plastic ear-tips (medium & large)
- Comfort-sealing ear-tips
- Multiple aeromedical certifications for helicopter winching / evacuation operations

DIMENSIONS

- Pkg: H 1.75 in. x W 6.5 in. x D 13 in.
- Device Length: L 30.5 in.
- Weight: 13.5 oz



Learn more about this and other Products with a Mission® at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277