

WEARABLE NECK WOUND SIMULATOR

REV071123

Supports Active Bleeding Simulation with Live Actors

The Wearable Neck Wound Simulator from North American Rescue Simulation is a silicone prosthetic device enabling live action role players to simulate a real-world neck injury. The soft, flexible, light toned skin wraps comfortably around the neck and attaches with the sewn hook and loop connectors. This realistic injury features exposed tissue and muscle and includes an integrated bleeding system tube that allows role players to demonstrate active bleeding during scenarios. Enhancing realism creates a valuable training event that better prepares first responder trainees for real-world situations. Easily cleaned and reusable, this tool enables rapid preparation for training exercises and minimizes moulage requirements – while maximizing realism.

This wearable simulator is compatible with all TOMManikin® and TAMikin® simulators.

Specifications:

- Hook and loop closure
- Exposed tissue and muscle
- Integrated bleeding system

Dimensions:

- L 6 in. x W 6 in. x H 6 in.
- Weight: 1 lbs.

Wearable Neck Wound Simulator

ITEM #	NSN #
93-0215	----



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**