

# IV POLE *Adjustable*

TACEVAC



## Features:

- Essential, lightweight litter accessory designed to fit the pre-molded holes of our Talon® and Raven® framed litters
- Easy to transport and fully adjustable, from a minimum clearance of 18 in. expanding to a height of 32 in. when deployed

## Efficient and Reliable Soution...

North American Rescue's Adjustable IV Pole is exclusively designed to expand the treatment capabilities of our family of evacuation platforms by fitting the pre-molded holes of all NAR framed litter platforms including the Talon®, Raven®, T3, Medevac & OSL litters.

This essential, lightweight accessory is easy to transport and fully adjustable, from a minimum clearance of 18 in. expanding to a height of 32 in. when deployed.

*Exclusively designed to fit the pre-molded holes of our framed litters and weighs less than 5 ounces.*

## Specifications:

- Retracted: H 19 in. x W 8.75 in. x D 0.6 in.
- Fully Extended Height: 33 in.
- Weight: 4.8 oz.

## Adjustable IV Pole

Item#	NSN#	CE
60-0008	6530-01-504-9054	

*An essential and lightweight accessory designed to expand treatment capabilities. Easy to transport and fully adjustable.*

PRODUCTS WITH A  
**MISSION**®



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

NOTE: All accessories are sold separately. Litters are not included.