

LITTER ACCESSORIES

Buckle Replacement Set or Retention Strap

Litter Support Bar Retention Strap...

Designed to prevent the support bar from collapsing when moving the Talon or Raven litters over obstacles. The strap is made of 1 in. military grade nylon with an ITW side release buckle with an adjustable tension feed.

Features:

- Prevents support bar from collapsing when moving
- Made of 1 inch Military Grade Nylon
- Adjustable Tension Feed

Specifications:

- Packaged:
L 3.6 in. x W 1.25 in. x D 1.6 in.
- Extended: L 104 in. x W 1 in.
- Weight: 2 oz



Retention Strap

Item#	NSN#
60-0026	



Litter Buckle Replacement Set

State of the art dual adjust model ITW Nexus TSR 200-BLK features patented Slick-Click latching mechanism and contoured profiles for maximum webbing grip while allowing smooth and easy adjustment. One package includes one male and one female component.

This buckle set is compatible with all NAR's litters (Talon, Raven, T3, MedEvac and OSL) as well as the NAR Rescue Sleeve II. Buckle set is only available in color Black.

Features:

- State of the art dual adjust model ITW Nexus TSR 200-BLK
- Packaged with one male and one female component
- Slick-Click latching mechanism
- Contoured profiles for maximum webbing grip
- Smooth and easy adjustment

Specifications:

- Packaged:
L 3.75 in. x W 2.375 in.
x D 0.5 in.
- Weight: 1 oz



Buckle Replacement Set

Item#	NSN#
ZZ-0032	6530-01-435-7337



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®