

Learn more about this and
other Products with a Mission®
at www.NARescue.com

FLAT EMERGENCY TRAUMA DRESSINGS - ETD®

FLAT FOLDED - 4 & 6 IN.

Flat folded rapid, controlled bandaging

The Flat Emergency Trauma Dressing (ETD®) series of bandages provide an effective, constant pressure bandage used for the treatment of wounds to control bleeding. The Flat ETD® is packaged in a low cube space, flat folded configuration to reduce its diameter and footprint in kits. It consists of a resilient elastic wrap equipped with a sterile non-adherent pad, durable securing device and innovative Quick-Grip Roll Control configuration. This combination facilitates a rapid, controlled, and extremely effective application process in the most demanding conditions.

The ETD® can also be used to sling and swath an arm, to secure splints to an extremity or to perform in any other function requiring an elastic wrap.

4 IN. FLAT ETD® DIMENSIONS

- Flat Package: L 4.25 in. x W 2.75 in. x D 0.875 in.
- Bandage: L 56 in. x W 4 in.
- Pad: H 6 in. x W 4 in.
- Weight: 2.6 oz



FEATURES/SPECIFICATIONS

- Suitable for applying immediate direct pressure
- Flat folded, vacuum configuration reduces cube space and footprint in your kit
- Multi-functional, all-in-one design
- Simple application that is easy to remember under combat/survival stress

- Omission of pressure bars and hooks to moderate the use of fine motor skills
- Quick-Grip Roll Control tabs to prevent the bandage from unrolling during application
- Vacuum-sealed, Low-Cube Packaging, with NAR's signature Red-Tip Technology® red tear notches



6 IN. FLAT ETD® DIMENSIONS

- Flat Pkg: L 6 in. x W 3.5 in x D 0.875 in.
- Bandage: L 56 in. x W 6 in.
- Pad: H 6 in. x W 6 in.
- Weight: 4 oz



FLAT EMERGENCY TRAUMA DRESSINGS - ETD® (FLAT FOLDED)

SIZE STYLE	ITEM #	NSN #
4 in. Flat	30-0093	---
6 in. Flat	30-0094	---



Flat Folded Configuration



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277