

EMERGENCY TRAUMA DRESSINGS - ETD®

STANDARD ROLLED - 4 & 6 IN.

Learn more about this and other Products with a Mission® at www.NARescue.com



Rapid, controlled, effective bandaging

The Emergency Trauma Dressing (ETD®) series of bandages are engineered using data from over 63 After Action Reviews (AARs) of real-world battle dressing applications. The ETD® consists of a resilient elastic wrap equipped with a sterile non-adherent pad, durable securing device and innovative Quick-Grip Roll Control configuration. This combination facilitates a rapid, controlled, and extremely effective application process in the most demanding conditions.

The ETD® can also be used to sling and swath an arm, to secure splints to an extremity or to perform in any other function requiring an elastic wrap.

4 IN. ETD® SPECIFICATIONS

- Standard Package: L 4.5 in. x W 2.5 in. x D 1.75 in.
- Bandage: L 56 in. x W 4 in.
- Pad: H 6 in. x W 4 in.
- Weight: 2.6 oz



6 IN. ETD® SPECIFICATIONS

- Standard Package: L 6 in. x W 2.5 in. x D 2.1 in.
- Bandage: L 56 in. x W 6 in.
- Pad: H 6 in. x W 6 in.
- Weight: 4 oz



FEATURES/SPECIFICATIONS

- Suitable for applying immediate direct pressure
- Multi-functional, all-in-one design
- Simple application that is easy to remember under combat/survival stress
- Omission of pressure bars and hooks to moderate the use of fine motor skills
- Quick-Grip Roll Control tabs to prevent the bandage from unrolling during application
- Vacuum-sealed, low-cube packaging, with NAR's signature Red-Tip Technology® red tear notches

EMERGENCY TRAUMA DRESSINGS - ETD® (STANDARD)

| SIZE STYLE | ITEM # | SEL # | NSN # |
|----------------|---------|--------------|------------------|
| 4 in. Rolled | 30-0031 | 09MS-04-BAND | 6510-01-558-4108 |
| 6 in. Rolled | 30-0032 | 09MS-04-BAND | 6510-01-558-4114 |



STERILE



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277