

EMERGENCY TRAUMA DRESSING - ETD®

STANDARD ROLLED - 4 & 6 IN.

REV072823

4 in. ETD® Specifications:

- Std. Pkg: L 4.5 in. x W 2.5 in. x D 1.75 in.
- Bandage: L 70 in. x W 4 in.
- Pad: H 6 in. x W 4 in.
- Weight: 2.6 oz



6 in. ETD® Specifications:

- Std. Pkg: L 6 in. x W 2.5 in. x D 2.1 in.
- Bandage: L 70 in. x W 6 in.
- Pad: H 6 in. x W 6 in.
- Weight: 4 oz



Rapid, controlled, effective bandaging

The Emergency Trauma Dressing (ETD®) series of bandages are engineered using data from over 63 After Action Reviews (AARs) of real-world battle dressing applications. The ETD® consists of a resilient elastic wrap equipped with a sterile non-adherent pad, durable securing device and innovative Quick-Grip Roll Control configuration. This combination facilitates a rapid, controlled, and extremely effective application process in the most demanding conditions.

The ETD® can also be used to sling and swath an arm, to secure splints to an extremity or to perform in any other function requiring an elastic wrap.

Features:

- Suitable for applying immediate direct pressure
- Multi-functional, all-in-one design
- Simple application that is easy to remember under combat/survival stress
- Omission of pressure bars and hooks to moderate the use of fine motor skills
- Quick-Grip Roll Control tabs to prevent the bandage from unrolling during application
- Vacuum-sealed, low-cube packaging, with NAR®'s signature Red-Tip Technology® red tear notches

Emergency Trauma Dressings - ETD® (Standard)



Size Style	Item#	SEL#	NSN#
4 in. Std.	30-0031	09MS-04-BAND	6510-01-558-4108
6 in. Std.	30-0032	09MS-04-BAND	6510-01-558-4114



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**