

NAR[®] MINI RESPONDER - ETD[®]

EMERGENCY TRAUMA DRESSING

Our most compact compression bandage designed to provide an effective pressure bandage when reduced packaging size is required...

The NAR[®] Mini Responder ETD[®] from North American Rescue[®] is a flat folded 42 in. elastic bandage with 4 x 6 inch non-adherent pad designed to be easily stored on the person in pant pockets, shirt pocket, trauma plate pouches, etc. It's small, compact packaging is less than an inch thick with easy open Red-Tip Technology[®] tear notches. The Mini Responder ETD[®] has both hook & loop as well as a sturdy C-Clip closure to hold the dressing in place. The bandage is sterile unless opened or damaged.

DIMENSIONS:

- Packaged:
L 5 in. x W 3.5 in. x D 0.75 in.
- Pad: 4 in. x 6 in.
- Bandage length (unstretched): 42 in.
- Weight: 2.4 oz



NAR[®] MINI RESPONDER ETD[®]

WIDTH	ITEM #	NSN #
4 in.	30-0179	---



STERILE

FEATURES/SPECIFICATIONS

- Small, compact packaging
- Flat folded
- Both hook and loop as well as C-Clip securing device
- 4 x 6 inch non-adherent sterile pad
- 42 inch elastic bandage (unstretched)
- Sterile unless package is opened or damaged



Learn more about this and other **Products with a Mission[®]** at www.NARescue.com



NORTH AMERICAN RESCUE[®]
www.NARescue.com • 888.689.6277