

Individual Kit (CCRK)

Designed with all the Tools for the Warfighter/Operator to Mitigate Extremity Bleeding, Tension Pneumothorax and Airway Obstruction.

The Individual Kit from North American Rescue® has been designed for the Warfighter/Operator and packed to provide all of the tools to mitigate the three leading causes of preventable combat death (60% extremity bleeding / 33% tension pneumothorax / 6% airway obstruction).

Our modular system can be attached to the operator using the split-tale connector, PALS/MOLLE style webbing or by fastening to the web/ gun belt using the quick weave method. Belt and leg attachment straps can be removed and reattached as needed to allow maximum flexibility during load carriage configuration.



FEATURES/SPECIFICATIONS

- 500D CORDURA® Solution-Dyed Nylon with Near Infrared Signature Reduction
- Ghillie-Tex IR signature reduction hardware
- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- MOLLE Style Attachment System Compatible
- Removable Belt and Leg Straps
- Leg Strap with Quick-Grip™
- Available in NAR®'s 11 Point Berry Amendment Nylon Standard



DIMENSIONS

- Case Closed: H 8.5 in. x W 8.75 in. x D 4.5 in.
- Weight: 1 lb 12 oz

KIT CONTENTS

- 1 x Bag (CCRK® - Individual)
- 4 x Black Talon® Nitrile Trauma Gloves (2 pr.)
- 1 x NAR® Nasopharyngeal Airway with Lube 28F (117mm)
- 1 x HyFin® Chest Seal
- 1 x ARS® Needle Decompression Kit (14 G x 3.25 in.)
- 1 x C-A-T® (Combat Application Tourniquet®) - Black
- 1 x ETD® 6 in. Emergency Trauma Dressing
- 2 x NAR® S-Rolled Gauze (4.5 in. x 4.1 yd)
- 1 x Polycarbonate Eye Shield
- 1 x Combat Casualty Reference Card
- 1 x Triage Card (NAR T2 Tag®)

INDIVIDUAL KIT (CCRK) R_X ONLY

ITEM #	COLOR	OPTION
80-0003	Black	6530-01-529-1134
80-0004	OD Green	---
80-0005	Coyote	---



Learn more about this and other Products with a Mission® at www.NARescue.com