

# SQUAD KIT

*Enhance the survivability of the assault squad*

Developed with the end-user requirements and skill set in mind, our innovative Combat Casualty Response Kit® - Squad Kit™ capitalizes on durability, versatility, mobility, and efficiency.

The kit's medical supplies are carefully selected to address the top three causes of death from penetrating trauma as identified by the CoTCCC (60% extremity bleeding / 33% tension pneumothorax / 6% airway obstruction).

Multiple-load carriage configuration capabilities include fanny pack, over the shoulder, attachments for 3-day assault packs and MOLLE style webbing.

## FEATURES/SPECIFICATIONS

- 500D CORDURA® Solution-Dyed Nylon with Near Infrared Signature Reduction
- Ghillie-Tex IR signature reduction hardware
- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- Multiple-load carriage configuration capabilities: fanny pack, over the shoulder, attachments for 3-day assault packs and MOLLE style webbing
- Enhanced Hip Stabilization System for greater security during transport
- Trauma Shear Storage on Front for Quick Access

## SQUAD KIT (CCRK)<sup>Rx ONLY</sup>

ITEM #	COLOR	OPTION
80-0037	Black	6545-01-543-5721
80-0038	OD Green	6545-01-543-3159
80-0039	Coyote	6545-01-543-2956
80-0251	MTC	---

## AVAILABLE COLORS



## DIMENSIONS

- Case Closed: H 7.75 in. x W 13 in. x D 9 in.
- Weight: 4 lb 6.5 oz

*Contoured profile, large central compartment, multiple attachment points & front trauma-shear storage*



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

# SQUAD KIT

Enhance the survivability of the assault squad



## KIT CONTENTS

- 1 x Bag (CCRK® Squad)
- 5 pr Black Talon® Nitrile Trauma Gloves
- 2 x NAR® Nasopharyngeal Airway with Lube - 28F (117mm)
- 1 x Petrolatum Gauze (3 in. x 18 in.)
- 1 x HyFin® Chest Seal
- 2 x ARS® Needle Decompression Kit (14G x 3.25 in.)
- 2 x C-A-T® (Combat Application Tourniquet®) - Black
- 4 x ETD® 6 in. Emergency Trauma Dressing
- 4 x NAR® S-Rolled Gauze (4.5 in. x 4.1 yd)
- 1 x ETD® Abdominal Emergency Trauma Dressing
- 2 x SAM® Splint II
- 1 x NAR® Trauma Shears (7.25 in.)
- 4 x Polycarbonate Eye Shield
- 1 x Surgical Tape (2 in.)
- 1 x Combat Casualty Reference Card
- 2 x Triage Card (NAR T2 Tag®)



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)