

MARINER KIT

Designed to Meet the Rigorous Requirements of the Emergency Med Assist Team (EMAT)

North American Rescue's Mariner Kit is the standard casualty response bag for the Submarine Forces.

The compact bag is orange with a reflective trim for high visibility, easy transport and convenient storage. The large central compartment, multiple attachment points and front Trauma-shear storage offer a versatile, mobile platform with ready-to-use supplies for treating leading causes of preventable combat death.

The Mariner Kit is packed in a contoured profile bag with a large central compartment, multiple attachment points and front trauma shear storage. The kit's contoured profile is designed to be worn as a belly bag or slung over the shoulder.

FEATURES/SPECIFICATIONS

- Constructed of durable 1000 Denier Nylon with reflective trim
- Sand & water resistant zipper with silent pulls
- Multiple-load carriage configuration capabilities: fanny pack, over the shoulder, attachments for 3-day assault packs and MOLLE style webbing
- Trauma Shear Storage on front for quick access



Standard Casualty Response Bag for the Submarine Forces

DIMENSIONS

- Case Closed: H 7.75 in. x W 13 in. x D 9 in.
- Weight: 6 lbs



KIT CONTENTS

- 1 x Bag (Mariner)
- 5 pr Black Talon® Nitrile Gloves - Black, Lg
- 2 x NAR® Nasopharyngeal Airway with Lube - 28F (117mm)
- 1 x Petrolatum Gauze (3 in. x 18 in.)
- 1 x HyFin® Vent Chest Seal-Twin Pk
- 2 x C-A-T® - Black (Combat Application Tourniquet®) - Black
- 4 x ETD® 6 in. Emergency Trauma Dressing
- 4 x NAR® S-Rolled Gauze (4.5 in. x 4.1 yd.)
- 1 x ETD® Abdominal Emergency Trauma Dressing
- 2 x SAM® Splint II
- 2 x NAR® Saline Lock Kit
- 1 x BOA® IV Constricting Band
- 1 x NAR® Trauma Shears (7.25 in.)
- 4 x Polycarbonate Eye Shield
- 1 x Surgical Tape (2 in.)
- 1 x Combat Casualty Reference Card
- 2 x Triage Card (NAR T2 Tag®)

MARINER KIT Rx ONLY

ITEM #	COLOR	OPTION
80-0107	ORG	6545-01-543-5770

Learn more about this and other Products with a Mission® at www.NARescue.com

