

BAC-Pack™ Bougie-Aided Cricothyroidotomy Pack

A comprehensive solution that allows the user to perform a surgical cric in any environment

The North American Rescue® Bougie-Aided Cricothyroidotomy Pack - BAC-Pack™ is a comprehensive solution that allows the user to perform a surgical cric in any environment. With the components arranged by order of use and organized with elastic loops, the BAC-Pack is an easy-to-use kit in a small pouch made of durable 500D nylon.

The BAC-Pack™ utilizes a purpose-built tube that has a pre-inserted bougie for rapid application, and the bougie is fully adjustable based on the patient's or provider's needs. The excess bougie can also be conveniently stored in the one-of-a-kind bougie cap/stopper. The tube meets the TCCC recommendations for the open surgical technique (flanged and cuffed airway cannula of less than 10 mm outer diameter, 6-7 mm internal diameter, and 5- 8 cm of intratracheal length).

The BAC-Pack™ also comes with the latest innovation in safety scalpels, the #10 Protected Scalpel, for effortless penetration through the cricothyroid membrane and the NAR® Tracheal Hook to maintain the opening without worry. Two alcohol pads and a gauze pad are included for skin prep, and a 10cc syringe and tube securing strap are included for finishing the procedure.

FEATURES/SPECIFICATIONS

- All-in-one cricothyroidotomy pack for use in any environment
- Pre-inserted bougie in tube allows for a more rapid, easy procedure
- Tube meets TCCC recommendations for bougie-aided open surgical technique
- #10 Protected Scalpel included for safety and ease of use
- Small, durable pouch keeps the items organized and ready to use

DIMENSIONS

- Packaged: H 8 in. x W 4 in. x D 1.75 in.
- Weight: 5.2 oz

BAC-PACK™

ITEM #	NSN
10-0052	6515-01-732-3942

SINGLE USE ONLY 
 RX ONLY  
NOT MADE WITH NATURAL RUBBER LATEX



Learn more about this and other Products with a Mission® at www.NARescue.com




BAC-Pack™ Bougie-Aided Cricothyroidotomy Pack



Kit Contents:

- 1x Rolled Nylon Pouch
- 2x Alcohol Prep Pad
- 1x Gauze Pad (2 pack) - 4 in. x 4 in.
- 1x #10 Protected Scalpel
- 1x NAR® Tracheal Hook
- 1x NAR® Tracheostomy Kit with Bougie Introducer (Cuffed, 6mm)
- 1x Tracheal Tube Securing Strap
- 1x 10cc IV Syringe



Featuring the advantage of NAR®'s Tracheostomy Kit with Bougie Introducer for assistance in proper tube placement

Learn more about this and other Products with a Mission® at www.NARescue.com

