

WHEELED BASE *for PAXLight®*

MED/SURGE

Made so PAXLight can be Independently Supported and Wheeled into Position...

Your PAXLight assembly can be independently supported and wheeled into appropriate position. No matter the evacuation platform or hospital bed utilized this allows you an easy independent support for your PAXLight.

Wheeled Base equipped with low center of gravity weight to secure your PAXLight during use. Also, comes with convenient attachment point for PAXLight power cord and AC adapter. Includes threaded self locating pin to secure PAXLight to Wheeled Base.

Wheeled base utilizes a 38 in. wide, 5 caster wheel base design and counter weight for stability

Features:

- Wheeled base utilizes a 38 in. wide, 5 caster wheel base design with 13.5 lb counter weight to stabilize the PAXLight assembly
- Also, included is a 29 in. long, 3/4 inch OD threaded rod that attaches directly to the weight block

Dimensions:

- Base: W 38 in.
- Assembled: H 38 in.
- Weight:
Base - 13.5 lb
Assembled - 22 lb



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

Wheeled Base for PAXLight®

Item#	NSN#
60-0045	

Threaded self locating pin secures the PAXLight to the wheeled base.

PRODUCTS WITH A
MISSION®