

EMS RAPID DEPLOYMENT KIT

REV092023



Effective treatment of casualties in unconventional environments...

The EMS Rapid Deployment Kit is designed to give pre-hospital healthcare providers the highest quality trauma supplies for the effective treatment of casualties in unconventional environments.

This kit and its contents are compact, durable and allow reliable and versatile interventions wherever the need arises. The container is constructed of the finest materials coupled with stringent sewing standards for maximum performance under the most extreme conditions.

The kit can be carried around the waist as a belly-bag or slung over the shoulder in bandolier fashion. The level of careful detail in its design and configuration sustains the science and experience gained from the valuable lessons learned in the streets and on the urban battlefield.

With its large central compartment, aggressive elastic storage system and multiple attachment points, the rugged EMS Rapid Deployment bag is also available for individual purchase.

Features:

- Gives pre-hospital healthcare providers mobility at a mass casualty event
- Incorporates lessons learned from urban battlefields
- Designed for efficient use under survival stress
- Premium container construction of durable 1000D CORDURA® nylon with reflective trim

EMS Rapid Deployment Kit

COLOR	ITEM #	NSN #
Orange	80-0151	----



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

EMS RAPID DEPLOYMENT KIT

REV092023



This kit and its contents are compact, durable and allow reliable and versatile interventions wherever and whenever the need arises.

Dimensions:

- Case Closed: H 7.75 in. x W 13 in. x D 9 in.
- Weight: 4 lb 12.8 oz

Kit Contents:

- 2 x C-A-T® (Combat Application Tourniquet®) - Orange
- 4 x NAR® S-Rolled Gauze
- 4 x NAR® Emergency Trauma Dressing, 6 in.
- 1 x NAR® Abdominal Emergency Trauma Dressing
- 2 x NAR® Pre-Lubricated Nasopharyngeal Airway, 28F
- 1 x Cyclone® New & Improved BVM
- 1 x Tactical Suction Device
- 1 x HyFin® Vent Chest Seal - Twin Pack
- 2 x ARS® for Needle Decompression, 14G
- 2 x SAM® Splint II
- 2 x Polycarbonate Eye Shield
- 1 x NAR® Surgical Tape (2 in.)
- 5 pr Bear Claw® Nitrile Gloves - Sand, Lg
- 1 x NAR® Trauma Shears, Lg
- 2 x Triage Card (NAR® T2 Tag®)



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**