

FIELD CORPSMAN KITS *BASIC OR DELUXE*

REV092123



PICTURED:
Deluxe Kit

Lightweight and Durable ENT Diagnostic Sets for Field Care...

The Basic Field Corpsman Kit from NAR® is a lightweight, durable and portable diagnostic kit provides excellent performance and flexibility in the field.

This water resistant and highly durable kit includes a variety of high intensity lights and magnifiers for day or night examinations.

Choose from the Basic or Deluxe portable diagnostic kit.
(Deluxe Kit contains all the items from the Basic kit plus a Mini Ophthalmoscope, additional Spare Bulb and extra Batteries).



Field Corpsman Kits

OPTION	ITEM#	NSN#
Basic	20-0001	6515-01-458-6057
Deluxe	20-0002	6515-01-458-6178



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**



Basic & Deluxe Kit Contents:

- General examination light: LED lamp to provide a focused beam of bright white light for superior illumination. Lightweight, water resistant and guaranteed unbreakable (batteries included)
- Red lens: Fits on the exam light for night use
- Fiber optic light bender/nasal head: To effectively direct light into the nasal or oral cavity. Can be used to transilluminate the sinuses
- 2 spare AAA batteries
- Soft case with mesh pouches: Secures the exam lights and accessories. Developed with padded sides and durable zipper. Two additional mesh pouches are positioned on the inside compartment for more parts or supplies
- Color instruction booklet: Contains operating instructions and part numbers
















Deluxe Kit Contents

includes everything from the Basic Kit above PLUS:

- Heine Mini 3000 Ophthalmoscope (batteries and spare bulb included)

Dimensions:

- BASIC:
L 8 in. x W 4.25 in. x D 2.25 in.
Weight: 6.2 oz
- DELUXE:
L 8 in. x W 4.25 in. x D 2.25 in.
Weight: 11.6 oz

 QTY 1 Pocket Exam Light	 QTY 1 Nasal Head Light Bender	 QTY 1 Otoscopic Head	 QTY 1 Speculum Adapter
 QTY 10 Disposable Adult Specula	 QTY 1 Eye Magnet	 QTY 1 Nylon Eye Loop	 QTY 1 Lighted Stainless Curette
 QTY 1 Cobalt Blue Lens	 QTY 1 Red Lens	 QTY 2 Fluor-I Strips	 QTY 2 Spare AAA Batteries
Deluxe Also Includes			
 QTY 1 Spare Ophthal. Bulb	 QTY 2 Spare AA Batteries	 QTY 1 Heine Mini Ophthalmoscope	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**



Offers healthcare professionals in the field excellence and flexibility as they care for the eyes and ears of patients.

Ear Care:

- Otoscope head: Provides bright illumination and 3X magnification for excellent visualization. Adjustable head allows examiners to adjust for personal comfort
- Reusable adult specula: For use with the otoscopic head
- Disposable adult specula: For use with the otoscopic head. Oscope accepts standard disposable specula
- Lighted stainless curette: Used to remove wax from the ear or to assist in debris removal from small wounds

Eye Care:

- Cobalt blue lens: Provides UV illumination to examine the eyes for corneal abrasions or small foreign bodies. Use with or without fluorescent eye dyes
- High resolution eye loop: Soft nylon loop fits on the cobalt blue lens. Designed to gently sweep foreign bodies from the surface of the eye. The lighted tip illuminates work area
- Lighted eye magnet: For the removal of metal foreign bodies in the eye
- Fluor-I-Strips: For use with cobalt blue lens for identification of corneal abrasions or to locate small particles in the eye



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**