

# RAPTOR™ MULTI-TOOL by Leatherman®

REV072524

## A Must Have Multi-Tool...

North American Rescue® is proud to offer the latest innovation from Leatherman®, the industry leader in application-driven, quality-made multi-tools for the last 30 years.

With the right balance of features for medical-specific ops without an overkill of options, the Leatherman® Raptor™ is a high-performance multi-tool equipped with a variety of mission critical tools essential to emergency medical practitioners in the field.

Components include precision 420HC stainless steel medical shears with a stow-away strap cutter and built-in oxygen tank wrench that may be quickly accessed with one hand, and can be used when the shears are open or closed. Additional components include an ultra-strong ring cutter that will easily cut through most standard gold bands, a carbide-tip glass breaker and handy 5 cm ruler.

The foldable shears come apart for easy cleaning and feature serrated blades to help grip materials for a strong fast cut. The ergonomic handles also feature glass-filled nylon scales and sport a lanyard hole and replaceable pocket clip for a variety of useful attachment points.

The Raptor™ comes complete with an injection molded polymer holster which is also MOLLE-compatible and allows device storage in the open or folded position. The holster's secure rotating clip easily attaches to your vest, pack or duty belt allowing rapid access to your multi-tool when you need it.

Backed by a 25 year warranty, the Leatherman® Raptor™ is built to last... a must have for every mobile medic and 1st responder working in life-threatening emergency situations.

### Leatherman® Raptor Multi-Tool

COLOR	ITEM #	NSN #
■ Black	70-0014	5110-01-626-8766
■ BLU/BLK	70-0264	----
■ RED/BLK	70-0266	----
■ BLU/YLW	70-0326	----
■ ORG/BLK	70-0327	----



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**

WITH A **MISSION®**

# RAPTOR MULTI-TOOL by Leatherman®

REV072524

Disassembles for  
easy cleaning.



## Features:

- **420HC Stainless Steel Medical Shears 1** with serrations to help grip materials for a strong fast cut
- One-handed **Strap Cutter 2** with safety lock liner works whether the shears are open or closed
- Stow-away **Oxygen Tank Wrench 3** stays out of the way until needed
- Ultra-strong **Ring Cutter 4** works on most standard gold bands
- Carbide-tip **Glass Breaker 5** quickly breaks standard civilian-grade glass during stressful rescue situations
- Convenient embedded 5 cm **Ruler 6**
- Built-in attachment points: **Lanyard Hole 7** & replaceable **Pocket Clip 8**
- Device folds down for small cube space
- Shears disassemble for easy cleaning & maintenance
- Includes injection-molded polymer **Holster 9** with rotating clip for rapid access from a vest, pack or duty belt; allows device storage in the open or folded positions; also MOLLE-compatible

## Dimensions:

- Shear length: 1.9 in.
- Closed length: 5 in.
- Weight: 5.8 oz

 **LEATHERMAN®**  
Leave nothing undone.™



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**