

NAR BP/STETHOSCOPE Combo Kit

MED/SURGE



Ultra-Compact Kit can be Stored in Most Trauma Kits...

This ultra-compact, multi-pocketed kit contains a set of 4 blood pressure cuffs and stethoscope combo which are instantly accessible and can be stored in most trauma kits.

This combo kit features a Stethoscope from ADC with extra large bell and ultra-sensitive diaphragm coupled with a set of four rugged nylon cuffs. The cuff set features adult, large adult, child and infant size cuffs and one-tube bladder with quick release luer slip connector, as well as an ambidextrous palm gauge with high impact casing and trigger valve.

*Kit contains ultra-compact blood-pressure cuffs
Set of 4 (adult, large adult, child, infant) and ADC stethoscope.*



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

NAR BP/Stethoscope Combo Kit

Item#	NSN#
20-0057	6515-01-542-7723

Includes a multi-pocketed and heavy duty zippered black nylon case, labeled for easy cuff size identification.

PRODUCTS WITH A
MISSION®

Specifications:

- SEL#: 09ME-03-STET
- Ultra-compact blood-pressure cuffs
- Set of 4 (adult, large adult, child, infant)
- One-tube bladder with quick-release luer slip-connector
- Ambidextrous palm gauge with high-impact casing and trigger valve
- ADC stethoscope with extra-large bell and ultra-sensitive diaphragm

Kit Contents:

- Ultra-compact blood pressure cuffs in a set of 4 (adult, large adult, child, infant)
- Ambidextrous palm gauge with high-impact casing and trigger valve
- One-tube bladder with quick-release luer slip-connector
- ADC stethoscope with extra-large bell and ultra-sensitive diaphragm
- Multi-pocketed, heavy-duty, zippered nylon case

Dimensions:

- Case Closed - L 6.75 in. x W 9 in. x H 5 in.
- Weight: 2 lb



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS WITH A
MISSION®