

# FIELD BLOOD TRANSFUSION KITS FBTK

REV092823

***“The advantage of the direct transfusion of human blood in cases of severe haemorrhage which we encounter in the emergencies of military surgery cannot to our minds, be overestimated.”<sup>1</sup>***

*The above comment, from a paper presented by Lieutenant Colonel A. Primrose of the Canadian army Medical Corps in 1916, remains as pertinent today as during those earlier days of military transfusion medicine...*

The Field Blood Transfusion Kit (FBTK) from North American Rescue® is a comprehensive kit providing all the needed supplies to allow the user to rapidly perform a field blood transfusion. The Transfusion of fresh whole blood is the best option for the management of casualties with hemorrhagic shock. Blood transfusions have been used successfully in military conflicts and other austere conditions to resuscitate trauma patients suffering from shock due to hemorrhage. Developed with extensive research and feedback from the Armed Services Blood Program Program Office the Fresh Whole Blood Transfusion Kit contains the supplies and materiel needed to test for ABO/Rh compatibility and perform a fresh whole blood transfusion in a sealed package.

<sup>1</sup>Primrose A, Ryerson ES. The Direct Transfusion of Blood: Its Value in Haemorrhage and Shock in the Treatment of the Wounded in War. Br Med J. 1916;2:384–386.

## **Specifications:**

- Comprehensive kit allowing the user to rapidly perform a field blood transfusion
- All supplies are compatible with the Joint Theater Trauma System Clinical Practice Guidelines
- Equipped with Donor/Recipient Modules to verify blood types of the donor and recipient
- Also allows collecting of blood samples from the donor for additional surveillance testing

## **Dimensions:**

- L 20 in. x W 13 in.
- Weight: 1.42 lbs.



## **Field Blood Transfusion Kits** Rx ONLY

ITEM #	NSN #
83-0000	6515-01-707-2935



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
 WITH A **MISSION®**

# FIELD BLOOD TRANSFUSION KITS *FBTK*

REV092823

## Donor Module Kit Contents:



### Donor Module Kit:

- 1 x Surgical Tape, 1 in.
- 2 x Alcohol Prep Pad
- 2 x Gauze Pad, 4 in. x 4 in. - Pk of 2
- 1 x Permanent Marker - Black, Lg
- 1 pr Bear Claw® Nitrile Gloves - Sand, Lg
- 1 x Donor Placard
- 2 x Vacutainer - Red Top, 6 ml
- 1 x CDP-A Single Blood Bag, 450 ml
- 1 x Kelly Curved Hemostat Disposable, 5.5 in.
- 1 x ASBP 572-EWB V29 Card
- 1 x IV Constricting Band - White, 1 in. x 18 in.
- 4 x Vacutainer - Purple Top, 6 ml
- 1 x Clear Adh. Dressing, 4 in. x 5 in.
- 1 x 21 ga Vacutainer Holder - Eclipse w/ Needle, 1-1/4 in.
- 1 x Blood Bag Sticker
- 1 x 550 Cord, 10 in.
- 1 x Beaded Security Cord, 11 in. (Marked at 6.5")
- 2 x Swabsticks - Pk of 3

## Recipient Module Kit Contents:



### Recipient Module Kit:

- 1 x Surgical Tape, 1 in.
- 2 x Alcohol Prep Pad
- 2 x Gauze Pad, 4 in. x 4 in. - Pk of 2
- 1 x Permanent Marker - Black, Lg
- 1 x Eldon Card
- 1 pr Bear Claw® Nitrile Gloves - Sand, Lg
- 1 x Recipient Placard
- 1 x IV Set Y-Type w/ Filter 170-260 Um w/o Pump
- 1 x SF-518 Card
- 1 x IV Constricting Band - White, 1 in. x 18 in.
- 1 x Clear Adh. Dressing, 4 in. x 5 in.
- 1 x 16 ga IV Catheter, 1.25 in.
- 1 x 18 ga IV Catheter - Safety, 1.25 in.
- 2 x Swabsticks - Pk of 3



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
 WITH A **MISSION®**