

MEDIC TRAUMA SHEET BAG ONLY

REV100423



Features:

- Versatile, rugged, lightweight equipment access platform
- Clear windows for easy viewing
- Oversized pockets
- Large grommet attachment points
- Quick attach-and-release straps secure to the skin of an aircraft or vehicle
- Bag only: Contents and carabiners not included

Specifications:

- Bag Open: H 41.5 in. x W 21.5 in. x D 2.5 in.
- Bag Folded: H 20.25 in. x W 21.5 in. x D 5 in.
- Weight: 4 lb 12 oz



Easy Access...

Easy access to life-saving equipment is what North American Rescue's Trauma Sheet is all about. **"Easy access"** means this versatile, rugged, lightweight equipment access platform can streamline the management of combat casualties during CASEVAC operations.

New and improved design 2 piece construction allows for higher flexibility in small aircrafts and confined spaces.

"Easy access" includes visibility: Our gear bag has clear windows.

"Easy access" includes stability: The Trauma Sheet has quick-attach-and-release straps allowing it to be secured to the skin of an aircraft or vehicle until it is needed.

"Easy access" includes convenience: Our Trauma Sheet has rugged, oversized pockets and large grommet attachment points for quick setup and equipment retrieval under stressful conditions. Plus, the Trauma Sheet comes empty so contents can be customized to a medic's specifications.

"Easy access" includes size: The bag opens to a usable 41.5 by 21.5 by 5 inches and weighs only 4.5 pounds.

"Easy access" is what medics want, so "easy access" is what North American Rescue gives them.

MEDIC TRAUMA SHEET

COLOR	ITEM #	NSN#
BLK	80-0067	6515-01-540-7612
ODG	80-0068	6515-01-540-7585
COY	80-0069	6515-01-540-7603

BLK

ODG

COY



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®