

MULTI MISSION EXPEDITIONARY RESPONSE KIT MMERK

REV111423

The standard Multi-Mission Expeditionary Response Kit from North American Rescue® which meets the requirements for the Expeditionary Force Senior Medics is the customizable modular solution ideal for personalized gear organization spanning multiple Expeditionary Mission profiles. Designed to meet the rigorous requirements of the U. S. Navy Expeditionary Force Senior Medics. This rugged 3200 cubic inch backpack style bag features 6 detachable clear vinyl internal pouches in various sizes which can be self-labeled for the specific mission sets. It also has multiple internal MOLLE straps covered with hook-and-loop and versatile Multi-mission belt all to facilitate mission specific configurations tailored to suit individual user preferences. Additionally, the Multi-mission Trauma Pack converts to a convenient trauma sheet ideal for use in a vehicle, boat, or aircraft. Its dedicated loops provide the capability of rappelling or lowering the bag to its desired location.

Specifications:

- Customizable modular solution ideal for personalized gear organization spanning Multiple Mission Expeditionary profile
- 3200-cubic inch backpack with (6) detachable clear vinyl internal pouches in various sizes
- Includes Multi-mission belt for different configurations
- Converts to trauma sheet for vehicles, boats, or aircraft
- Internal features include MOLLE straps covered in hook and loop and channeling for electrical communications wiring
- Dedicated loops, providing rappelling or lowering of the bag

Dimensions:

- H 25.5 in. x W 11 in. x D 9 in.
- Weight: 25 lbs. 12 oz.

MMERK

Item #	NSN#
85-2420	----



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277


PRODUCTS
WITH A **MISSION®**

MULTI MISSION EXPEDITIONARY RESPONSE KIT *MMERK*

REV111423

MMERK Contents:

- 3 x C-A-T® - Black (Combat Application Tourniquet®)
- 6 x NAR® Compressed Gauze
- 4 x Combat Gauze®, Z-Fold Hemostatic Dressing
- 5 x NAR® ETD®, 4 in. (Emergency Trauma Dressing®)
- 2 x NAR® ETD® ABD/Stump, 6 in.
- 3 x NAR® Pre-lubricated Nasopharyngeal Airway, 28F
- 3 x HyFin® Vent Chest Seal - Twin Pack
- 2 x ARS® for Needle Decompression, 14G
- 2 x SPEAR for Needle Decompression, 10G
- 1 x Tactical Suction Device Hand Held
- 1 x Cyclone® New & Improved BVM
- 1 x I-Gel Supraglottic Airway, Size 4
- 2 x Saline IV Solution, 500cc
- 2 x Baxter IV Tubing
- 1 x Blood Transfusion Field Kit
- 1 x Intraosseous EZ-IO Manual Needle
- 1 x Blizzard Survival Blanket
- 2 x NAR® Survival Blanket
- 1 x Slishman Traction Splint - Compact
- 3 x SAM® Splint II
- 4 x Triangular Bandage - ODG 45 in. x 45 in. x 63 in.
- 4 x Burntec® Burn Dressing - 4 in. x 4 in.
- 2 x Burntec® Burn Dressing - 8 in. x 8 in.
- 1 x Burntec® Burn Dressing - 8 in. x 16 in.
- 10 x Gauze Cover Sponge, 4 in. x 4 in.
- 10 x Antiseptic Towelette
- 10 x Flexible Fabric Bandages
- 2 x Surgical Tape Roll, 1 in.
- 2 x NAR® PPE Kit™ (Personal Protection Equipment Kit)
- 1 x Microshield CPR Mask
- 1 x Hand Sanitizer
- 6 pr Bear Claw® Gloves - Sand, Lg
- 1 x Junctional Emergency Treatment Tool™ (JETT®)
- 4 x Polycarbonate Eye Shield
- 1 x Stethoscope (Littman)
- 1 x MightySat Pulse Oximeter
- 1 x Digital Thermometer
- 4 x AA Batteries
- 1 x First light TORQ with batteries
- 1 x Sharps Container
- 2 x NAR® Mini Duct Tape, 2 in. x 100 in.
- 2 x Self Grip Tape, 2 in.
- 1 x NAR® Trauma Shears, Lg
- 1 x Scissor Leash™
- 1 x Model 7 Rescue Hook
- 1 x Combat Casualty Reference Card
- 5 x Combat Casualty Response Card
- 2 x Permanent Marker - Black, Lg
- 1 x Medical Storage Case, Armadillo - Black



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**