

I-ROK KIT *Individual Rugged Outdoor Kit*

REV062524

Essential First-Aid Protected In A Personal Utility Case...

Every good adventure requires that everyone gets home safely. The new I-ROK (Individual Rugged Outdoor Kit) from North American Rescue® is the perfect companion for any outdoor environment. The I-ROK contains essential first-aid supplies to treat common injuries to include basic injuries such as cuts, scrapes, bites & stings as well as the life-saving equipment to control major bleeding.

The I-ROK Kit comes in the ultimate personal utility case that is crushproof, dustproof and adventure proof! Your first-aid and trauma supplies will be protected from the outdoor elements to help ensure your kit is ready if ever needed. Waterproof up to a depth of 2 meters for up to 30 minutes, IP68 rated protection from water, dirt, snow & dust and rubberized protective bumpers to protect against drops. The case handles can be used to secure the kit for your outing as you can strap it to your gear, ATV/UTV, motorcycle, camper or whatever you use to get you to where the adventures are. The hard case comes in a variety of colors and keeps your contents protected and secure.

Features:

- Contains essential first-aid supplies to treat basic injuries as well as control major bleeding
- R40 storage case is crushproof, dustproof and waterproof up to a depth of 2 meters
- Flexible lid organizer helps keep smaller medical supplies neatly stored
- Built-in pressure valve keeps water and dust out while balancing air pressure
- I-ROK hard case comes in a variety of colors and with the option of Basic or with Bleeding Control Dressing (hemostatic bandage) versions

Dimensions:

- Measurements (Interior): H 1.9 in. x W 7.6 in. x D 4.7 in.
- Measurements (Exterior): H 2.9 in. x W 9.8 in. x D 6.1 in.
- Weight: 2.41 lb



I-ROK Individual Rugged Outdoor Kit

COLOR	LEVEL	ITEM #
Black	Basic	82-0031
OD Green	Basic	82-0033
Orange	Basic	82-0034
Black	with BCD	82-0035
OD Green	with BCD	82-0037
Orange	with BCD	82-0038
Black	with Celox™	82-0066
OD Green	with Celox™	82-0068
Orange	with Celox™	82-0069






















NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH A **MISSION®**

I-ROK KIT *Individual Rugged Outdoor Kit*

REV062524

BASIC -or-
BASIC with Hemostatics:

<p>BASIC</p>  <p>QTY 1 Responder Compressed Gauze</p>	-OR-	<p>BASIC with Hemostatics</p>  <p>QTY 1 QuikClot® BCD</p>	-OR-	<p>BASIC with CELOX™ Gauze</p>  <p>QTY 1 5 CELOX™ Hemostatic</p>
 <p>QTY 1 C-A-T® Tourniquet</p>	 <p>QTY 1 NAR® Mini Responder 4" ETD®</p>	 <p>QTY 1 HyFin® Vent Compact, Twin Pk</p>	 <p>QTY 5 Butterfly Bandage</p>	 <p>QTY 5 1" x 3" Bandage</p>
 <p>QTY 3 Safetec Burn Gel</p>	 <p>QTY 1 Microshield CPR Mask</p>	 <p>QTY 3 Antibiotic Ointment</p>	 <p>QTY 2 Sting & Bite Antiseptic</p>	 <p>QTY 1 Survival Blanket</p>
 <p>QTY 3 Hand Sanitizer</p>	 <p>QTY 3 Salty Britches® Chafing Ointment</p>	 <p>PAIR 1 Responder Blue Gloves, LG</p>	 <p>QTY 1 NAR® Mini Duct Tape, 2 in.</p>	
		 <p>QTY 1 Permanent Marker, Large</p>	 <p>QTY 1 Instructions Bleeding Control</p>	



Kit Contents :

- 1 x C-A-T® Tourniquet - Orange
- 1 x NAR® Responder Compressed Gauze
- OR-
- QuikClot® BCD
- OR-
- CELOX™ Gauze
- 1 x NAR® Mini ETD®, 4 in. (Emergency Trauma Dressing)
- 1 x HyFin® Vent Compact Chest Seal - Twin Pk
- 5 x Butterfly Bandage
- 5 x Bandage, 1 in. x 3 in.
- 3 x Burn Gel
- 1 x MicroShield CPR Face Mask
- 3 x Antibiotic Ointment
- 2 x Sting & Bite Packet
- 1 x NAR® Survival Blanket
- 3 x Hand Sanitizer
- 3 x Salty Britches® Chafing Ointment, Single Use
- 1 pr Responder Nitrile Gloves, Large
- 1 x NAR® Mini Duct Tape, 2 in.
- 1 x Permanent Marker - Black, Large
- 1 x Bleeding Control Instruction Card
- 1 x Pelican™ Ruck R40 Case



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH A **MISSION®**