

I-ROK XL *Individual Rugged Outdoor Kit*

Essential First-Aid Protected In A Personal Utility Case...

Every good adventure requires that everyone gets home safely. The new I-ROK XL (Individual Rugged Outdoor Kit) from North American Rescue is the perfect companion for any outdoor environment. The I-ROK XL contains the same essential first-aid supplies of the Standard I-ROK Kit with extra supplies to accommodate multiple injuries. These kits are designed to treat common injuries including basic injuries such as cuts, scrapes, bites & stings as well as the life-saving equipment to control major bleeding.

The I-ROK XL comes in the larger ultimate personal R60 utility case that is crushproof, dustproof and adventure proof! Your first-aid and trauma supplies will be protected from the outdoor elements to help ensure your kit is ready if ever needed. Waterproof up to a depth of 2 meters for up to 30 minutes, IP68 rated protection from water, dirt, snow & dust and rubberized protective bumpers to protect against drops. The case handles can be used to secure the kit for your outing as you can strap it to your gear, ATV/UTV, motorcycle, camper or whatever you use to get you to where the adventures are. This ample hard case comes in a variety of colors and keeps your contents protected and secure.

Features:

- Contains essential first-aid supplies to treat basic injuries as well as control major bleeding with extra supplies for multiple injuries
- I-ROK XL features a larger R60 storage case that is crushproof, dustproof and waterproof up to a depth of 2 meters
- Flexible lid organizer helps keep smaller medical supplies neatly stored
- Built-in pressure valve keeps water and dust out while balancing air pressure
- I-ROK XL larger hard case comes in a variety of colors and with the option of Basic, with Bleeding Control Dressing, or with CELOX™ Hemostatic Dressing



REV121223

Dimensions:

- H 8.9 in. x W 5.4 in. x D 3 in.
- Weight: 3 lb.

I-ROK XL Individual Rugged Outdoor Kit

COLOR	LEVEL	ITEM #
Black	Basic	82-0048
OD Green	Basic	82-0050
Orange	Basic	82-0051
Black	with Bleeding Control Dressing	82-0052
OD Green	with Bleeding Control Dressing	82-0054
Orange	with Bleeding Control Dressing	82-0055
Black	with Celox™ Hemostatic Dressing	82-0070
OD Green	with Celox™ Hemostatic Dressing	82-0072
Orange	with Celox™ Hemostatic Dressing	82-0073



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

I-ROK XL Individual Rugged Outdoor Kit


REV121223




QTY 1 C-A-T® Tourniquet



BASIC QTY 1 Responder Compressed Gauze



-OR- BASIC with Bleeding Control Dressing QTY 1 QuikClot® BCD



-OR- BASIC with CELOX™ Hemostatic Dressing QTY 1 5" CELOX™ Hemostatic



QTY 1 NAR® Mini Responder 4" ETD®



QTY 1 HyFin® Vent Compact, Twin Pk



QTY 2 4" x 4" Gauze Pad 2 PK



QTY 10 Butterfly Bandage



QTY 10 1" x 3" Bandage



QTY 1 BurnTec® 4"x4" Dressing



QTY 1 Triangular Bandage



QTY 5 Safetec Burn Gel



QTY 1 Microshield CPR Mask



QTY 5 Antibiotic Ointment



QTY 3 Sting & Bite Antiseptic



QTY 3 Salty Britches® Chafing Ointment



QTY 5 Hand Sanitizer



PAIR 2 Responder Blue Gloves, LG



QTY 1 Survival Blanket



QTY 1 Splinter Forceps, 4.5"



QTY 1 NAR Trauma Shears, Lg



QTY 1 NAR® Mini Duct Tape, 2 in.



QTY 1 Permanent Marker, Large



QTY 1 Instructions for Bleeding Control



Kit Contents:

- 1 x C-A-T® Tourniquet, Orange
- 1 x NAR® Responder Compressed Gauze
- OR-
- 1 x QuikClot® Z-Fold Bleeding Control
- OR-
- 1 x CELOX™ First Aid Hemostatic
- 1 x NAR® Mini Responder ETD®, 4in. (Emergency Trauma Dressing)
- 1 x HyFin® Vent Compact Chest Seal - Twin Pack
- 10 x Butterfly Bandage
- 10 x Bandage, 1 x 3 in.
- 2 x Gauze Pad, 4 x 4 in. (2 Pk)
- 1 x BurnTec® Dressing, 4 x 4 in.
- 5 x Burn Gel
- 1 x MicroShield CPR Face Mask
- 5 x Antibiotic Ointment
- 3 x Sting & Bite Packet
- 1 x NAR® Survival Blanket
- 5 x Hand Sanitizer
- 3 x Salty Britches® Chafing Ointment - Sgl
- 2 pr Responder Nitrile Gloves - Blue, Large
- 1 x NAR® Mini Duct Tape, 2 in.
- 1 x Permanent Marker - Black, Large
- 1 x Bleeding Control Instruction Card
- 1 x Triangular Bandage
- 1 x Splinter Forceps, 4.5 in.
- 1 x NAR® Trauma Shears, Large
- 1 x Pelican Ruck R60 Case

The same essential first-aid supplies from our I-ROK Kit with EXTRA supplies to accommodate multiple injuries...



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**