

MARITIME ASSAULT KIT

Complete With All the Tools to Help Mitigate the Three Leading Causes of Preventable Combat Death

North American Rescue's Maritime Assault Kit™ meets the unique and specific requirements of the Navy Expeditionary Combat Command for a small, efficient IFAK.

Complete with all the tools to mitigate the three leading causes of preventable combat death (60% extremity bleeding / 33% tension pneumothorax / 6% airway obstruction).

The Maritime Assault Kit Modular system can be attached to the operator in three ways:

1. Using the split-tale connector with height adjustment and Quick-Grip™ leg strap
2. PALS/MOLLE style webbing for attachment to the vest/body armor
3. Fasten to the web/gun belt using the quick weave method. Belt and leg attachment straps can be removed and reattached as needed to allow maximum flexibility during load carriage configuration

FEATURES/SPECIFICATIONS

- 500D CORDURA® Solution-Dyed Nylon with Near Infrared Signature Reduction
- Ghillie-Tex IR signature reduction hardware
- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- MOLLE-style Attachment System Compatible
- Removable Belt and Leg Straps
- Leg Strap with Quick-Grip™
- Available in NAR®'s 11 Point Berry Amendment Nylon Standard



DIMENSIONS

- Case Closed: H 7.4 in. x W 6.9 in. x D 5 in.
- Weight: 1 lb. 13 oz.

MARITIME ASSAULT KIT Rx ONLY

ITEM #	COLOR	OPTION
80-0013	Black	6545-01-529-4185
80-0014	OD Green	6545-01-531-6626
80-0015	Coyote	6545-01-529-4187
80-0245	MTC	---

AVAILABLE COLORS



Learn more about this and other Products with a Mission® at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

MARITIME ASSAULT KIT



KIT CONTENTS

- 1x Bag (CCRK® - Maritime)
- 1x C-A-T® (Combat Application Tourniquet®) - Black
- 1x ETD® 6 in. Emergency Trauma Dressing
- 1x Combat Gauze® Z-Fold Hemostatic Dressing
- 1x NAR® S-Rolled Gauze (4.5 in. x 4.1 yd)
- 1x NAR® Nasopharyngeal Airway with Lube - 28F (117mm)
- 1x HyFin® Vent Chest Seal, Twin Pack
- 1x ARS® for Needle Decompression (14G x 3.25 in.)
- 1x Polycarbonate Eye Shield
- 2 pr Black Talon® Nitrile Trauma Gloves
- 1x Combat Casualty Reference Card
- 1x Triage Card (NAR T2 Tag®)



Learn more about this and other Products with a Mission® at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277