

# MARITIME ASSAULT KIT

## Complete With All the Tools to Help Mitigate the Three Leading Causes of Preventable Combat Death

North American Rescue's Maritime Assault Kit™ meets the unique and specific requirements of the Navy Expeditionary Combat Command for a small, efficient IFAK.

Complete with all the tools to mitigate the three leading causes of preventable combat death (60% extremity bleeding / 33% tension pneumothorax / 6% airway obstruction).

 QTY 1 C-A-T® Tourniquet	 QTY 1 ETD® 6 in. Dressing	 QTY 1 NAR® Gauze S-Rolled	 QTY 1 Combat Gauze® Z-Fold
 QTY 1 Pre-lubed NAR® Naso, 28F	 QTY 1 HyFin® Vent Chest Seal, Twin Pk	 QTY 1 14G ARS® Needle Decomp.	 PAIR 2 Black Talon Nitrile Gloves, LG
 QTY 1 PES Eye Shield	 QTY 2 7 Hook Safety Cutter	 QTY 1 T2® Tag Triage Card	 QTY 1 Casualty Reference Card



## The Maritime Assault Kit Modular system can be attached to the operator in three ways:

1. Using the split-tale connector with height adjustment and Quick-Grip™ leg strap
2. PALS/MOLLE style webbing for attachment to the vest/body armor
3. Fasten to the web/gun belt using the quick weave method. Belt and leg attachment straps can be removed and reattached as needed to allow maximum flexibility during load carriage configuration

## Specifications:

- Case Closed: H 7.4 in. x W 6.9 in. x D 5 in.
- Weight: 1 lb. 13 oz.

## MARITIME ASSAULT KIT (CCRK)

COLOR	ITEM#	NSN#
BLK	80-0013	6545-01-529-4185
ODG	80-0014	6545-01-531-6626
COY	80-0015	6545-01-529-4187
MTC	80-0245	

Rx ONLY

## Available Colors:



**NORTH AMERICAN RESCUE**<sup>®</sup>  
www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A **MISSION**<sup>®</sup>

# MARITIME ASSAULT KIT

REV120123

## Features:

- 500D CORDURA® Solution-Dyed Nylon with Near Infrared Signature Reduction
- Ghillie-Tex IR signature reduction hardware
- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- MOLLE-style Attachment System Compatible
- Removable Belt and Leg Straps
- Leg Strap with Quick-Grip™
- Available in NAR®'s 11 Point Berry Amendment Nylon Standard

## Kit Contents:

- 1 x Bag (CCRK®-Maritime)
- 1 x C-A-T® (Combat Application Tourniquet®)
- 1 x ETD® 6 in. Emergency Trauma Dressing
- 1 x Combat Gauze® Z-Fold Hemostatic Dressing
- 1 x NAR® S-Rolled Gauze (4.5 in. x 4.1 yd)
- 1 x NAR® Pre-Lubricated Nasopharyngeal Airway, 28F
- 1 x HyFin® Vent Chest Seal, Twin Pack
- 1 x ARS® for Needle Decompression (14G x 3.25 in.)
- 1 x Polycarbonate Eye Shield
- 2 pr Black Talon® Nitrile Trauma Gloves
- 1 x Combat Casualty Reference Card
- 1 x Triage Card (NAR T2 Tag®)



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**