

# PAX COOLER BAG

## FOR THE CRÉDO PROMED SERIES FOUR 4 LITER CONTAINER

Introducing the North American Rescue® PAX Cooler Bag for the Crêdo ProMed Series Four 4 Liter Container. Carry both your temperature controlled medical supply needs as well as your blood & IV fluid warming device in the same bag for rapid access.

The NAR® PAX Cooler Bag features a main compartment custom sized to fit the Thermal Insulation Chamber - TIC™ and VIP components from the Crêdo ProMed Series Four 4 Liter Container. Simply transfer the removable TIC™ and VIP's into the NAR® PAX Blood Cooler Bag from your Crêdo.

In addition to the temperature control compartment, the bag features two detachable pouches that are ideal for storing your blood and IV fluid warming device and other essential IV access equipment.

Use the quick release buckles to detach the pouches from the bag. The NAR® PAX Bag is easy to transport using the reinforced nylon handles or the built-in shoulder strap.

The NAR® PAX Bag features the non-porous PAX Tech material that provides an easy-to-clean and decontaminate exterior when compared to traditional woven nylon materials.

Be prepared to easily transport Blood and IV fluids for field deployment with the Blood Cooler Bag for the Crêdo ProMed Series Four 4 Liter Container.

### FEATURES/SPECIFICATIONS

- Features the non-porous PAX Tech material that provides an easy-to-clean and decontaminate exterior
- Has two detachable pouches that are ideal for storing your blood and IV fluid warming device and other essential IV access equipment
- Compatible with the Thermal Insulation Chamber - TIC™ and VIP components from the 4 Liter Crêdo cooler

### PAX COOLER BAG FOR THE CRÉDO PROMED SERIES FOUR 4 LITER CONTAINER

ITEM #	NSN
35-0033	----

### DIMENSIONS

- H 8in. x W 11 in. x D 16 in.
- Weight: 2 lbs. 8 oz.



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277