

TRUECLOT® Tourniquet Trainer Simulator

REV081221

The TrueClot® Tourniquet Application Trainer is a wearable device that allows full application of a tourniquet for training purposes. Training blood can be pumped through the device to simulate hemorrhage, and when a tourniquet is correctly applied, bleeding is stopped. The padded arm cuff reduces pinching and pain and allows for proper application of a tourniquet during classroom or scenario training. The trainer can be worn under a uniform or large t-shirt.

Features:

- Durable ballistic nylon construction
- Easy to clean and air dry
- Apply tourniquet high and tight or 2" - 3" above the wound
- Heavy duty padding to reduce pinching and pain
- Bleeding wound with simulated vessel

Dimensions:

- H 12 in. x W 16 in. x D 7 in. (packaged)
- Weight: 5 lbs

TrueClot® Tourniquet Trainer Simulator

ITEM #	NSN#
93-0931	----

Proper tourniquet placement and application are crucial to effective hemorrhage control.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

TRUECLOT® Tourniquet Trainer Simulator

REV081221

Kit Contents:

- Trainer Sling Assembly (Left Arm Fit Only)
- Liter Squeeze Bottle
- Quick Disconnect Fitting
- Plastic Case



Ideal for training tourniquet placement techniques featuring a durable silicone construction with no moving parts, reinforced tubing mount, and a more realistic wound pattern and outer anatomy...



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**