

MEDIC TRAUMA PACK

Our Most Comprehensive Trauma Management Platform for Tactical Healthcare Professionals in the Line of Duty

The Medic Trauma Pack is a light framed platform with silent, dirt and sand-resistant zippers that close over a large main compartment with modular, removable pull bags for all your trauma management essentials.

The adjustable shoulder and sternum straps make for a comfortable fit. Cinch straps are included for the sides to maintain load integrity during stress maneuvers. Grommets are provided at the base of the pack for drainage.

Includes 9 removable, color-coded pull bags and 2 removable side compartments to provide modular organization for all your medical equipment. With over 1,000 cubic inches of gear capacity this is our most comprehensive kit/pack available.

Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Over 1000 cubic inches of gear capacity
- Contoured, padded shoulder straps to distribute load over a larger area for comfortable carry
- Flap allows for Hydration Tube to extend out of pack
- Constructed of durable 500 Denier Nylon
- Sand & water resistant zipper with silent cord pulls
- External Cinch Straps to tighten pack and decrease load shift
- MOLLE Style Attachment System Compatible

MEDIC TRAUMA PACK (CCRK) Rx ONLY

ITEM #	COLOR	OPTION
80-0070	Black	6545-01-572-7251
80-0071	OD Green	6545-01-594-5617*
80-0072	Coyote	6545-01-588-1862

*NATO NSN#

DIMENSIONS


- Case Closed: H 19 in. x W 19 in. x D 10 in.
- Weight: 18 lbs. 1 oz.



AVAILABLE COLORS



MEDIC TRAUMA PACK



KIT CONTENTS

- 2 x C-A-T® (Combat Application Tourniquet®) - Black
- 4 x NAR® S-Rolled Gauze (4.5 in. x 4.1 yd)
- 6 x ETD® 6 in. -Emergency Trauma Dressing
- 2 x ETD® Abdominal -Emergency Trauma Dressing
- 2 x NAR® Nasopharyngeal Airway with Lube - 28F (117mm)
- 1 x I-Gel Supraglottic Airway Device
- 1 x Cyclone® New & Improved BVM
- 2 x ARS® Needle Decompression Kit (14G x 3.25 in.)
- 1 x HRS™ Heat Reflective Shell
- 4 x SAM® Splint II
- 1 x CT-6 Tactical Traction Splint
- 1 x Tactical Suction Device
- 2 x Petrolatum Gauze (3 in. x 18 in.)
- 6 x Polycarbonate Eye Shield
- 2 x Surgical Tape (2 in.)
- 4 x NAR® Saline Lock Kit
- 1 x NAR® Sharps Shuttle
- 1 x BOA® IV Constricting Band
- 1 x NAR® Trauma Shears (7.25 in.)
- 5 pr Black Talon® Nitrile Trauma Gloves
- 2 x NAR® PPE Kit™ -Personal Protection Equipment Kit
- 1 x Armadillo™ Medication Storage Case
- 1 x Combat Casualty Reference Card
- 4 x Triage Card (NAR T2 Tag®)
- 1 x NAR® Headlamp
- 1 x CCRK® Medic Trauma Bag

Light Framed Platform and Removable Pull Bags for all your Trauma Management Essentials



Learn more about this and other Products with a Mission® at www.NARescue.com