

MINI TPAK TRAIL PERSONAL AID KIT

REV062524

Introducing the smaller sibling of our Trail PAK (Personnel Aid Kit)

The Mini Trail Personal Aid Kit (Mini T-PAK). Still featuring life saving items like its bigger sibling but in a compact aLOKSAK® package. The aLOKSAK® bag is comprised of a 6 Mil patent protected new generation film with a Hermetic Seal that ensures waterproofing and airtight seal. The bag is certified waterproof up to 200 Feet / 60 Meters. The Mini T-PAK equips you to address the first aid needs of injuries ranging from major bleeding trauma and fractures to minor cuts and abrasions. Be prepared to save a life by bringing the Mini T-PAK on all your adventures.

Ideal For:

- Backpacking
- Hiking
- Hunting
- Camping
- Vehicles
- Adventure travel
- Surfing
- Motocross / ATV
- Mountain Biking
- Skiing
- Snowboarding
- Individual Sports

Mini Trail Personal Aid Kit

LEVEL	ITEM #	NSN#
Basic	82-0075	----
with BCD	82-0076	----

Contents

QTY 1 C-A-T® Tourniquet

QTY 1 NAR® Gauze, Wound Packing

QTY 1 QuikClot® BCD

QTY 1 4" Flat Elastic Wrap

QTY 1 Microshield CPR Mask

QTY 1 HyFin® Vent Compact, Twin Pk

QTY 1 Survival Blanket

QTY 1 Minor Care Bandages Kit

QTY 3 Antiseptic Towelette

QTY 3 Antibiotic Ointment

QTY 3 Sting & Bite Antiseptic

QTY 1 IvyX Wipe Post Contact

QTY 2 Salty Britches® Chafing Ointment

QTY 1 Moleskin Bandage (6 Pk)

QTY 1 Gecko Grip® Flat Tape, 2"

PAIR 1 Responder Blue Gloves, LG

Features:

- Life saving equipment in an ultra-portable, waterproof, and re-sealable aLokSak® bag
- Addresses first aid needs of injuries ranging from major bleeding trauma to minor cuts and abrasions



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH A **MISSION®**