

H.I.D.E. SUSPENDERS *Hidden Internal Drag Extraction*

RESCUE

Everyday equipment support system with the added benefit of an integrated drag system...

North American Rescue's H.I.D.E. Suspenders (Hidden Internal Drag Extraction) rapidly deploys a convenient drag system that allows one or two rescuers to quickly move an injured person safely behind cover.

It's drag loop remains securely stored inside the long channels, but for training can be rapidly re-packed into the training loops. The system is constructed of MIL-W-17337 Mil-Spec nylon webbing which easily resets during training without having to be removed.

Features:

- Everyday equipment support system with the added benefit of an integrated drag system
- Rapidly deploying drag system allows one or two rescuers to quickly move an injured person behind cover
- Drag loop remains securely stored inside the long channels, but for training purposes can be rapidly re-packed into the training loops
- Constructed of MIL-W-17337 Mil-Spec nylon webbing

Dimensions:

- Drag Handle Deployed: 38 in.
- Size: One size fits most
- Weight: 11 oz
- Hardware: Minimum 200 lb straight pull
- Webbing: 1 in. (800 lb) | 1.5 in. (1,200 lb)



NOTE:
Suspenders only.
All accessories sold separately



H.I.D.E. Suspenders

Color	Item#	NSN#
BLK	90-0079	
ODG	90-0081	

Rapidly deploying drag system with allows a quick move of an injured person behind cover

PRODUCTS WITH A
MISSION[®]



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277