


# RIG SERIES EAGLE IFAK BAG

REV050224



## **Low profile IFAK bag offers rapid access to your C-A-T® tourniquet...**

The RIG Series Eagle IFAK is a new low profile individual first aid platform designed with a protective hood and quick latch release to keep your life-saving equipment protected and rapidly accessible for deployment. Highlighted design advancements include easy-one hand quick release hood latch and zippers, detachable inner shock cord panel that detaches from the IFAK bag and attachment system that allows both vertical (MOLLE backing) and horizontal (3 in. belt loop) carriage orientations.

### **Features:**

- Nylon bag with sewn in embroidered medical cross identifier
- Protective hood with quick latch release helps protect contents (not included)
- Allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Quick pull latch and tab for easy one-handed opening
- Detachable (hook & loop backed) shock-cord inner panel adjusts to secure contents (not included) & is easily removable for rapid access

### **Specifications:**

- Packaged: H 8 in. x W 4 in. x D 3 in.
- Weight: 7 oz

### **RIG Series Eagle IFAK**

Color	Item#	NSN#
Black	80-0434T	----
Coyote	80-0435T	----
Multi-Cam	80-0436T	----
Ranger Green	80-0988T	----



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277