

# LARYNGOSCOPE SET

## SPECIAL OPERATIONS

### Perform direct visual laryngoscopy even in the most severe environments

The lightweight, compact, and highly innovative Special Operations Laryngoscope Set from North American Rescue® is designed to operate effectively under the most extreme environmental conditions. Professionals rely on it. To ensure the laryngoscopes perform as long and as well as necessary, the sets are packaged with spare LED bulbs and batteries.

The fiber optic handle is water-proof, sand-proof, and dust-proof, and features a textured gripping surface. For convenience, the set includes single-use fiber optic blades made from an enhanced polymer for exceptional strength and rigidity. For added convenience, the handle's removable light source is designed to work with any standard (green-spec) fiber optic blade.

The complete set rolls up for easy storage and easy transport.

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



### FEATURES/SPECIFICATIONS

- Lightweight & compact
- Resistant to dust, sand and water
- Operates in severe environmental conditions
- Textured, water-proof & dust-proof handle
- 80 Lumen LED examination lamp
- Removable LED light source works with any standard fiber optic blade
- Single-use fiber optic blades of enhanced polymer for exceptional strength and rigidity

### SPECIAL OPS LARYNGOSCOPE SET

ITEM #	NSN
10-0001	6515-01-450-9790



### DIMENSIONS

- Open: L 8 in. x W 5.5 in. x D 2 in.
- Closed: L 8 in. x W 2 in. x D 2 in.
- Weight: 10 oz.



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

# LARYNGOSCOPE SET

## SPECIAL OPERATIONS



## KIT CONTENTS

- 1 x Laryngoscope Pocket Light™ with two AAA Batteries (Item# 10-0074)
- 1 x Laryngoscope Handle (Item# 10-0005 | NSN# 6515-01-460-4736)
- 1 x Miller #3 Laryngoscope Blade (Item# 10-0011 | NSN# 6515-01-460-4681)
- 1 x Mac #3 Laryngoscope Blade (Item# 10-0013 | NSN# 6515-01-460-4685)
- 1 x Laryngoscope Pouch (Item# 10-0006 | NSN# 6515-01-467-7624)
- 2 x spare AAA batteries
- 1 x Laryngoscope Set Booklet IFU (rev. 02/15/24)

*\*NOTE: To order replacement parts, please call Customer Service at 1-888-689-6277*

## SPECIAL OPS LARYNGOSCOPE SET

ITEM #	NSN
10-0001	6515-01-450-9790



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277