

NAR® TRAUMA SHEARS

Learn more about this and other Products with a Mission® at www.NARescue.com



Premium-grade stainless steel with precision sharpness...

Built for safety and efficiency, the NAR® Trauma Shears from North American Rescue let tactical healthcare professionals get the job done, and get it done right.

These paramedic shears are made from premium-grade stainless steel. Their serrated lower edge cuts through the toughest materials. Their soft-grip, non-slip contoured handles provide maximum control, even when wet. Their clothing lift lip and rounded tip provide the ultimate in patient safety.

NAR Trauma Shears are made to meet the challenges of casualty care with precision and dependability, their unique all-over black color makes them the shears of choice in tactical environments.

Our medical scissors come in two convenient lengths, 6.25 or 7.25 inches.

NAR TRAUMA SHEARS

ITEM #	OPTION	NSN
ZZ-0063	7.25 in.	6515-01-538-9276
ZZ-0064	6.25 in.	6510-01-622-1624*

* NATO NSN#

FEATURES/SPECIFICATIONS

- Premium-grade stainless steel with precision sharpness
- Enhanced serrated lower edge cuts through the toughest materials
- Clothing-lift lip and rounded tip for ultimate patient safety
- Soft-grip, non-slip contoured handles for max control even when wet
- Unique all-over black color
- Autoclave Safe 143°C / 290°F

Unique all-over black Trauma Shears ideal for the tactical environment



6.25 in. Small

7.25 in. Large featuring O₂ wrench

DIMENSIONS

6.25 in. Shears:

- H 6.25 in. x W 2.75 in. x D 0.3 in.
- Weight: 1.5 oz.

7.25 in. Shears:

- H 7.25 in. x W 3.7 in. x D 0.4 in.
- Weight: 2 oz.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277