

CHITOGAUZE® XR PRO



Introducing ChitoGauze® XR PRO:

A revolutionary hemostatic dressing for the external, temporary control of severely bleeding wounds. It is the next generation of chitosan dressing that is composed of a polyester/rayon blend non-woven medical gauze that is coated with chitosan. Featuring an embedded x-ray detectable strip, it is Z-folded and packaged in a vacuum sealed pouch that is compact and rugged.

ChitoGauze® XR PRO maximizes hemostatic performance and reduces the risk of re-bleeding by providing effective hemostasis outside of the normal clotting cascade. ChitoGauze® XR PRO also has natural antibacterial properties against a wide range of gram positive and gram negative microorganisms which are beneficial in presence of debris and environment before the injured are treated with antibiotics or are transferred to a medical treatment facility. Designed for the battlefield, ChitoGauze® XR PRO is an efficacious and easy-to-apply dressing that maximizes hemostatic performance before severe bleeding turns deadly.

Features:

- Z-Fold Vacuum Packed for easy application
- Featuring an embedded x-ray detectable strip
- Readily conforms into the wound surfaces with complex geometries to rapidly staunch bleeding and seal the wound site
- Offers antimicrobial properties against 26 bacterial organisms
- Non-woven, polyester/rayon blend medical gauze is coated with chitosan to maximize its performance and control severe bleeding

Specifications:

- Gauze: L 12 ft. x W 3 in.
- Pkg: H 6 in. x W 5 in. x D 0.65 in.
- Weight: 0.7 oz.

Item#	NSN#
30-0049	6510-01-591-7740



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS WITH A
MISSION®