

MINI MEDIC BAGS

Learn more about this and other Products with a Mission® at www.NARescue.com



Compact and Versatile for different mission profiles

The Mini Medic Bag is a compact trauma aid bag designed to provide various capabilities for different mission profiles. External pouches and side/bottom MOLLE allow for external tourniquet or other attachment options. Large pouch on back panel includes a drop-down pouch with abrasive resistant padded bottom allows storage of compact breaching tools (not included). External expansion zipper allows for increased load capacity if needed.

The inside of the bag has two small and one large, clear tear-away pouches for medical or other equipment supply storage. Along with two mesh inner pouches and an internal pouch on the inner lid, many storage options exist.

Breaching tool not included.



Features Drop-down pouch for Compact Breaching or Hydration Packs (not included)

FEATURES/SPECIFICATIONS

- Compact, multi-purpose trauma aid bag designed to provide versatile organization for different mission profiles
- External storage options include 2 pockets and extra drop-down pouch for breaching tools or hydration pack (not included)
- Expansion zipper increases load capacity
- MOLLE capability on the sides, bottom and inside back of the bag
- Internal pouch storage includes 3 clear tear-away, 2 main compartment mesh pouches and 1 larger inner lid pouch
- Bag only option

DIMENSIONS

- H 14 in. x W 11 in. x D 5 in. (Depth expandable to 7.25 in.)
- Weight: 3 lb

MINI MEDIC BAGS

ITEM #	COLOR	NSN #
80-0954T	Black	---
80-0966T	Coyote	---
80-955T	OD Green	6545-01-722-6568
80-0995T	Ranger Green	---
80-1230T	Multicam	---

* Explore our packed out Mini Medic kits at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277