

NAR® SURVIVAL BLANKET

SILVER

Balances Performance, Durability, and Size

The North American Rescue® Survival Blanket is the ideal solution for a small, lightweight survival wrap as it balances performance, durability and size. The wrap reflects radiant body heat and is constructed of a combination of mylar and tear-resistant polyethylene fabric material.

The blanket is completely wind and waterproof, further aiding in retaining warmth. The NAR® Survival Blanket is an essential piece of equipment in a compact package to help maintain body heat.

Learn more about this and other **Products with a Mission®** at www.NARescue.com



FEATURES/SPECIFICATIONS

- Lightweight, durable space blanket made of a combination of mylar and polyethylene, perfect for regulating body temperature
- Reflective silver color helps to put radiant heat back into the body
- Wind and waterproof, further aiding in retaining warmth
- 56 in. x 84 in. wrap for full body coverage

NAR® SURVIVAL BLANKET

ITEM #	NSN
ZZ-0038	4240-01-697-8869

DIMENSIONS

- Packaged: L 5 in. x W 3 in.
- Blanket: L 84 in. x W 56 in.
- Weight: 1.84 oz.

